



# National Institute of Standards & Technology

## Certificate of Analysis

### Standard Reference Material<sup>®</sup> 2975

#### Diesel Particulate Matter (Industrial Forklift)

This Standard Reference Material (SRM) is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic hydrocarbons (PAHs) and nitro-substituted PAHs (nitro-PAHs) in diesel particulate matter and similar matrices. In addition to certified, reference, and information values for selected PAHs and nitro-PAHs, reference or information values are provided for total extractable mass, particle-size distribution, and specific surface area; supplemental information on mutagenic activity is also provided. All of the chemical constituents for which certified, reference, and information values are provided are naturally present in the diesel particulate material. A unit of SRM 2975 consists of a bottle containing 1 g of the diesel particulate matter collected from an industrial diesel-powered forklift.

SRM 1975 Diesel Particulate Extract [1], which is a dichloromethane extract of the diesel particulate matter, was prepared from the same lot of material as SRM 2975. A second diesel particulate material, SRM 1650b Diesel Particulate Matter [2], which was originally issued in 1985, is representative of heavy-duty diesel engine particulate emissions.

**Certified Mass Fraction Values:** Certified values are provided for PAHs in Table 1 and nitro-PAHs in Table 2. A NIST certified value is a value for which NIST has the highest confidence in its accuracy, in that all known or suspected sources of bias have been investigated or taken into account [3]. The certified values for the PAHs and nitro-PAHs are based on the agreement of results obtained at NIST from two or more independent analytical methods [3].

**Reference Mass Fraction Values:** Reference values are provided for additional PAHs in Table 3 and nitro-PAHs in Table 4. Some PAHs and nitro-PAHs are listed more than once in the table depending on the extraction conditions used (see "Preparation and Analysis"). Reference values for total extractable mass and the particle-size distribution are provided in Tables 5 and 6, respectively. Reference values are noncertified values that represent the best estimates of the true value; however, the values do not meet the NIST criteria for certification and are provided with associated uncertainties that may reflect only measurement precision, may not include all sources of uncertainty, or may reflect a lack of sufficient statistical agreement among multiple analytical methods.

**Information Values:** Information mass fraction values are provided for additional PAHs of relative molecular mass 302 in Table 7 and specific surface area as determined by N<sub>2</sub> gas adsorption is provided in Table 8. An information value is considered to be a value that may be of use to the SRM user, but insufficient information is available to assess the uncertainty associated with the value or only a limited number of analyses were performed.

**Expiration of Certification:** The certification of **SRM 2975** is valid, within the measurement uncertainty specified, until **31 January 2022**, provided the SRM is handled and stored in accordance with the instructions given in this certificate (see "Instructions for Handling, Storage, and Use"). However, the certification is nullified if the SRM is damaged, contaminated, or otherwise modified.

**Maintenance of SRM Certification:** NIST will monitor this SRM over the period of its certification. If substantive technical changes occur that affect the certification before the expiration of this certificate, NIST will notify the purchaser. Registration (see attached sheet) will facilitate notification.

Carlos A. Gonzalez, Chief  
Chemical Sciences Division

Robert L. Watters, Jr, Director  
Office of Reference Materials

Gaithersburg, MD 20899  
Certificate Issue Date: 03 December 2013  
*Certificate Revision History on Last Page*

SRM 2975

Page 1 of 11

The coordination of the technical measurements leading to the certification of SRM 2975 was under the leadership of M.M. Schantz and S.A. Wise of the NIST Chemical Sciences Division.

Analytical measurements for the certification of SRM 2975 were performed by H.M. Bamford, D. Bezabeh, M. Lopez de Alda, D.L. Poster, L.C. Sander, M.M. Schantz, P. Schubert, and L. Walton of the NIST Chemical Sciences Division. Specific surface area, and porosity data were provided by P. Scheepers of the Department of Epidemiology at Katholieke Universiteit Nijmegen, Nijmegen, The Netherlands. The particle-size distribution data were provided by Honeywell, Inc., Clearwater, FL.

The diesel particulate material was provided by M.E. Wright of the Donaldson Company, Inc., Minneapolis, MN.

Statistical consultation was provided by N.A. Heckert, S.D. Leigh and M.G. Vangel of the NIST Statistical Engineering Division.

Support aspects involved in the issuance of this SRM were coordinated through the NIST Office of Reference Materials.

## INSTRUCTIONS FOR HANDLING, STORAGE, AND USE

**Handling:** This material is a naturally occurring diesel particulate material and contains constituents of known and unknown toxicities and mutagenicities; therefore, extreme caution and care should be exercised during its handling and use.

**Storage:** Store SRM 2975 in its original bottle at temperatures below 30 °C and keep away from direct sunlight.

**Use:** Prior to removal of subsamples for analysis, the contents of the bottle should be mixed thoroughly. The recommended minimum sample size is 100 mg.

## PREPARATION AND ANALYSIS<sup>(1)</sup>

**Sample Collection and Preparation:** The diesel particulate material used to prepare SRM 2975 was obtained from M.E. Wright of the Donaldson Company, Inc., Minneapolis, MN. The material was collected from a filtering system designed specifically for diesel-powered forklifts [4]. The diesel particulate material was received at NIST in a 55-gallon drum. The material was removed from the drum and homogenized in a V-blender for 1 h and then stored in polyethylene bags. A total of 13.7 kg of diesel particulate material was homogenized; a total of 5.65 kg of material was extracted for preparation of SRM 2975 [1] and the remaining diesel particulate material was bottled for distribution as SRM 2975.

### Polycyclic Aromatic Hydrocarbons (PAHs)

The general approach used for the value assignment of the PAHs in SRM 2975 consisted of combining results from analyses using various combinations of different extraction techniques and solvents, cleanup/isolation procedures, and chromatographic separation and detection techniques. This approach consisted of Soxhlet extraction and pressurized fluid extraction (PFE) using dichloromethane (DCM), toluene, or toluene/methanol mixture, cleanup of the extracts using solid-phase extraction (SPE), followed by analysis using the following techniques: (1) reversed-phase liquid chromatography with fluorescence detection (LC-FL) for analysis of isomeric PAH fractions isolated by normal-phase LC (i.e., multidimensional LC) and (2) gas chromatography/mass spectrometry (GC/MS) for analysis of the PAH fraction on four stationary phases of different selectivity (i.e., a 5 % [mole fraction] phenyl-substituted methylpolysiloxane phase, a 50 % phenyl-substituted methylpolysiloxane phase, a non-polar, extra low bleed proprietary phase, and a smectic liquid crystalline stationary phase).

Multiple sets of GC/MS results, designated as GC/MS (Ia and Ib), GC/MS (II), GC/MS (III), GC/MS (IVa, IVb, and IVc), and GC/MS (Va through Vd), were obtained using three columns with different selectivities for the separation of PAHs.

---

<sup>(1)</sup>Certain commercial instruments, materials, or processes are identified in this report to adequately specify the experimental procedure. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the instruments, materials, or processes identified are necessarily the best available for the purpose.

*GC/MS (Ia and Ib):* For GC/MS (I) analyses, duplicate subsamples of 100 mg from eight bottles of SRM 2975 were extracted with toluene:methanol (1:1 volume fraction) using PFE as described by Schantz et al. [5]. The excess volume in the PFE cells was filled with clean sodium sulfate. The extracts were concentrated to about 0.5 mL and placed on an aminopropylsilane SPE cartridge from which they were eluted with 20 mL of 2 % DCM in hexane (volume fraction). The eluant was concentrated and then analyzed by GC/MS using a 0.25 mm i.d.  $\times$  60 m fused silica capillary column with a 5 % (mole fraction) phenyl-substituted methylpolysiloxane phase (0.25  $\mu$ m film thickness; DB-5MS, Agilent Technologies, Wilmington, DE), and these results are designated as GC/MS (Ia). A subset of 8 of the 16 extracts from GC/MS (Ia) were also analyzed on a 50 % (mole fraction) phenyl-substituted methylpolysiloxane stationary phase (0.25 mm i.d.  $\times$  60 m, 0.25  $\mu$ m film thickness; DB-17MS, Agilent Technologies), and these results are designated as GC/MS (Ib).

*GC/MS (II through V):* For the GC/MS (II) analyses, 100 mg samples of SRM 2975 were extracted with DCM using PFE; the extracts were processed and analyzed as described above for GC/MS (Ia). GC/MS (III) was identical to GC/MS (II) except that Soxhlet extraction with DCM for 18 h was used instead of PFE. For the GC/MS (IV) analyses, subsamples of 40 mg to 100 mg of SRM 2975 were extracted with DCM using PFE and the extracts processed as described above for GC/MS (I). The processed extracts were then analyzed by GC/MS using three different columns: 5 % phenyl methylpolysiloxane [GC/MS (IVa)], 50 % phenyl methylpolysiloxane [GC/MS (IVb)], and a 0.2 mm i.d.  $\times$  25 m (0.15  $\mu$ m film thickness) smectic liquid crystalline phase [SB-Smectic, Dionex, Lee Scientific Division, Salt Lake City, UT; GC/MS (IVc)].

In addition, a study investigating the effect of increasing the temperature and pressure used for PFE on the extraction efficiency for PAHs was performed. In this study, the solvent used was toluene although a 9:1 toluene:methanol (volume fraction) was also evaluated. The PFE conditions used included: 100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa [methods GC/MS (Va) through GC/MS (Vd), respectively]. Following an SPE step similar to GC/MS (I) above, the processed extracts were analyzed using a non-polar, extra low bleed proprietary phase (0.25 mm i.d.  $\times$  60 m, 0.25  $\mu$ m film thickness).

*Reversed-Phase Liquid Chromatography with Fluorescence Detection (LC-FL):* For the LC-FL analyses, subsamples of approximately 200 mg from each of six bottles of SRM 2975 were Soxhlet-extracted for 20 h using 200 mL of DCM. The extracts were concentrated and processed through aminopropylsilane SPE cartridges as described above for the GC/MS analyses. The processed extract was further fractionated using normal-phase LC on a semi-preparative aminopropylsilane column ( $\mu$ Bondapak NH<sub>2</sub>, 9 mm i.d.  $\times$  30 cm, Waters Associates, Milford, MA) to isolate isomeric PAH fractions as described previously [6–8]. Four fractions were collected containing PAHs of relative molecular mass 178 and 202 (fraction 1), 228 (fraction 2), 252 and 276 (fraction 3), and 278 (fraction 4). All of the PAH fractions were analyzed using a 5  $\mu$ m particle-size polymeric octadecylsilane (C<sub>18</sub>) column (4.6 mm i.d.  $\times$  25 cm, Hypersil-PAH, Keystone Scientific, Inc., Bellefonte, PA) with wavelength-programmed fluorescence detection [7–9]. For all of the GC/MS and LC-FL measurements described above, selected perdeuterated PAHs were added to the diesel particulate samples prior to solvent extraction for use as internal standards for quantification purposes.

*PAH Isomers of Relative Molecular Mass ( $M_r$ ) 302:* For the determination of the relative molecular mass ( $M_r$ ) 302 isomers, samples of approximately 100 mg each were extracted using PFE at 100 °C with DCM as the extraction solvent. The extracts were then concentrated with a solvent change to hexane, passed through an aminopropyl SPE cartridge, and eluted with 40 mL of 10 % DCM in hexane (volume fraction). The processed extract was then analyzed by GC/MS using a 0.25 mm i.d.  $\times$  60 m fused silica capillary column with a 50 % phenyl-substituted methylpolysiloxane phase (0.25  $\mu$ m film thickness; DB-17MS). Perdeuterated dibenzo[*a,i*]pyrene was added to the diesel particulate matter prior to extraction for use as an internal standard.

*Homogeneity Assessment for PAHs:* The homogeneity of SRM 2975 was assessed by analyzing duplicate samples of 100 mg each from eight bottles selected by stratified random sampling. Samples were processed and analyzed as described above for GC/MS (I). Statistically significant differences among bottles were observed for the PAHs at the 100 mg sample size, and this source of uncertainty has been incorporated in the calculation of the uncertainty associated with the assigned values.

### **Nitro-Substituted Polycyclic Aromatic Hydrocarbons (Nitro-PAHs)**

SRM 2975 was analyzed at NIST for the determination of nitro-PAHs. Three sets of three, four, and five samples of SRM 2975 (~100 mg each) were spiked with the following perdeuterated nitro-PAHs for use as internal standards: 9-nitroanthracene-*d*<sub>9</sub>, 3-nitrofluoranthene-*d*<sub>9</sub>, 1-nitropyrene-*d*<sub>9</sub>, and 6-nitrochrysene-*d*<sub>11</sub>. The samples were extracted using PFE with DCM or toluene as the extraction solvent.

Following concentration, each sample was processed through an aminopropylsilane SPE cartridge using 40 mL of 20 % DCM in hexane (volume fraction). The concentrated eluant was then subjected to normal-phase LC using a semi-preparative amino/cyano-phase column with a mobile phase of 20 % DCM in hexane to isolate the nitro-PAH fraction. The nitro-PAH fraction was analyzed by GC with negative chemical ionization mass spectrometry (GC/NCI-MS) using a 0.25 mm i.d.  $\times$  30 m fused silica capillary column containing a 5 % diphenyl-substituted dimethylsiloxane phase (0.25  $\mu$ m film thickness) or a 50 % phenyl-substituted methylpolysiloxane phase (0.25  $\mu$ m film thickness). Methane was used as the NCI-MS collision gas. The effect of increasing the temperature and pressure used in PFE was also investigated at NIST for the nitro-PAHs using 100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa.

The concentration of 1-nitropyrene was also measured at the Department of Epidemiology at Katholieke Universiteit Nijmegen, Nijmegen, The Netherlands using GC/MS [10].

### **Value Assignment for PAHs and Nitro-PAHs**

The current value assignment of PAHs and nitro-PAHs in SRM 2975 is based on the measurements used when SRM 2975 was issued in 2000 and on recent additional measurements using different extraction conditions. Previous studies on solvent extraction of diesel particulate matter for determination of PAHs indicated that PFE at 100 °C removed higher quantities of higher relative molecular mass PAHs than does Soxhlet extraction [5]. Recent studies on the extraction of PAHs and nitro-PAHs from diesel particulate matter [11–13] have shown that using PFE at 200 °C removes higher quantities of some (not all) PAHs and nitro-PAHs than using PFE at 100 °C. As a result of these studies, value assignment for PAHs and nitro-PAHs in SRM 2975 is based on measurements using PFE at both 100 °C and 200 °C and on limited measurements using Soxhlet extraction. In cases where the quantities of the individual PAHs and nitro-PAHs determined do not change with the extraction method or temperature, the measurements were combined and the resulting values are denoted as certified values. When different results are obtained using PFE at 100 °C or at 200 °C, the values are reported for both temperatures, and they are denoted as reference values. These reference values should be considered as “method dependent” values because they are dependent on the extraction method and temperature.

### **Total Extractable Mass**

For the determination of total extractable mass, six subsamples of approximately 1 g to 2 g of SRM 2975 were Soxhlet-extracted for 18 h with DCM. The extract was concentrated to approximately 20 mL and then filtered to remove particulate matter. Aliquots of 100  $\mu$ L to 150  $\mu$ L were placed in tared aluminum foil pans; the DCM was evaporated until constant mass was obtained, and then the mass of the remaining residue was determined. The total extractable mass reference value is provided in Table 5.

### **Particle-Size Information**

Dry particle-size distribution measurements for SRM 2975 were obtained as part of a collaborative effort with Honeywell, Inc., Clearwater, FL. A Microtrac particle analyzer, which makes use of light-scattering techniques, was used to measure the particle-size distribution of SRM 2975 [14]. Briefly, a reference beam was used to penetrate a field of particles and the light that scatters in the forward direction from the field is measured. The particle size as a volume distribution is derived via a computer-assisted analysis. From these data, the total volume, average size, and a characteristic width of the particle-size distribution were calculated.

### **Specific Surface Area and Porosity**

The specific surface area and porosity were determined based on nitrogen gas adsorption measurements [15]. The gas adsorption measurements were performed on a NOVA-1200 instrument, Quantachrome Corp., Boynton Beach, FL at 77 K after the samples were outgassed for 24 h at 120 °C under vacuum. The nitrogen isotherms were analyzed using the Brunauer-Emmet-Teller (BET) equation [16] to obtain the surface area (Table 8) and the Barrett-Joyner-Halenda (BJH) method [17] to obtain the porosity. Based on the BJH method, SRM 2975 shows a wide distribution of mesopores, but with substantial outer area. The pore diameter of the particles in SRM 2975 range from 4 nm to 35 nm with the greater number of particles at about 20 nm.

### **Mutagenicity Activity**

The mutagenicity of solvent extracts of SRM 2975 has not been determined; however, the reference values for mutagenicity of SRM 1975 were determined by using the the *Salmonella typhimurium*/mammalian microsome mutagenicity assay and should be applicable to extracts of SRM 2975. SRM 1975 was prepared by Soxhlet extraction using DCM of multiple samples of SRM 2975. Reference values for mutagenic activity for SRM 1975 are reported in the Certificate of Analysis for SRM 1975 [1] and described in detail by Hughes et al. [18].

Table 1. Certified Mass Fraction Values for Selected PAHs in SRM 2975

	Mass Fraction <sup>(a)</sup> (mg/kg)
2-Methylphenanthrene <sup>(b,c,d,e)</sup>	2.22 ± 0.21
1-Methylphenanthrene <sup>(b,c,d,e,f,g)</sup>	0.923 ± 0.057
Triphenylene <sup>(e,h,i)</sup>	5.32 ± 0.24
Benzo[ <i>j</i> ]fluoranthene <sup>(g,i,j)</sup>	0.819 ± 0.093
Benzo[ <i>a</i> ]fluoranthene <sup>(e,g,h)</sup>	0.066 ± 0.007
Dibenz[ <i>a,h</i> + <i>a,c</i> ]anthracene <sup>(e,f,i)</sup>	0.523 ± 0.047
Picene <sup>(e,f,g,i)</sup>	0.902 ± 0.091
Dibenzo[ <i>a,e</i> ]pyrene <sup>(e,g)</sup>	0.599 ± 0.024
Dibenzo[ <i>b,k</i> ]fluoranthene <sup>(e,g)</sup>	2.54 ± 0.08

<sup>(a)</sup> The certified values are weighted means of the mass fractions from two to six analytical methods [19]. The uncertainty listed with each value is an expanded uncertainty about the mean [19,20], with coverage factor,  $k = 2$ , calculated by combining within method variances with a between method variance [21] following the ISO/JCGM Guide [22,23].

<sup>(b)</sup> GC/MS (Ia) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with toluene:methanol mixture (100 °C with 13.8 MPa)

<sup>(c)</sup> GC/MS (II) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(d)</sup> GC/MS (III) on 5 % phenyl-substituted methylpolysiloxane phase after Soxhlet extraction with DCM

<sup>(e)</sup> GC/MS (V) on non-polar extra low bleed proprietary phase after PFE with toluene at two temperatures and two pressures (100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa)

<sup>(f)</sup> GC/MS (IVa) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(g)</sup> GC/MS (IVb) on 50 % phenyl-substituted methylpolysiloxane phase of same extracts as GC/MS (IVa)

<sup>(h)</sup> GC/MS (IVc) on a smectic liquid crystalline phase of same extracts as GC/MS (IVa)

<sup>(i)</sup> LC FL of isomeric PAH fractions after Soxhlet extraction with DCM

<sup>(j)</sup> GC/MS (Ib) on 50 % phenyl-substituted methylpolysiloxane phase of selected extracts from GC/MS (Ia)

Table 2. Certified Mass Fraction Values for Selected Nitro-PAHs in SRM 2975

	Mass Fraction <sup>(a)</sup> (mg/kg)
9-Nitrophenanthrene <sup>(b)</sup>	0.466 ± 0.013
3-Nitrophenanthrene <sup>(c)</sup>	0.190 ± 0.007
2-Nitrofluoranthene <sup>(c)</sup>	0.231 ± 0.032
3-Nitrofluoranthene <sup>(b)</sup>	3.80 ± 0.24
1-Nitropyrene <sup>(b)</sup>	35.2 ± 2.2
7-Nitrobenz[ <i>a</i> ]anthracene <sup>(b)</sup>	3.57 ± 0.32
6-Nitrochrysene <sup>(b)</sup>	2.45 ± 0.33

<sup>(a)</sup> The certified values are weighted means of the mass fractions from two to six analytical methods [19]. The uncertainty listed with each value is an expanded uncertainty about the mean [19,20], with coverage factor,  $k = 2$ , calculated by combining within method variances with a between method variance [21] following the JCGM Guide [22,23].

<sup>(b)</sup> Three methods using GC/NCI-MS on either a 5 % phenyl-substituted or 50 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa) and four methods using GC/NCI-MS on a 50 % phenyl-substituted methylpolysiloxane phase after PFE with toluene at two temperatures and two pressures (100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa)

<sup>(c)</sup> Two methods using GC/NCI-MS on a 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa) and four methods using GC/NCI-MS on a 50 % phenyl-substituted methylpolysiloxane phase after PFE with toluene at two temperatures and two pressures (100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa)

Table 3. Reference Mass Fraction Values for Selected PAHs in SRM 2975  
Based on Extraction Method and Conditions

Extraction Conditions	Mass Fractions <sup>(a)</sup> (mg/kg)		
Soxhlet Extraction or PFE at 100 °C			
Naphthalene <sup>(b,c,d,e)</sup>	1.78	±	0.48
Biphenyl <sup>(e)</sup>	0.323	±	0.014 <sup>(f)</sup>
Fluorene <sup>(e,g,h,i)</sup>	0.440	±	0.071
Phenanthrene <sup>(b,c,d,e,g,h,i)</sup>	17.3	±	1.0
Anthracene <sup>(b,c,e)</sup>	0.039	±	0.003
3-Methylphenanthrene <sup>(c,d,e,g,h,i)</sup>	1.06	±	0.07
9-Methylphenanthrene <sup>(d,e,j)</sup>	0.448	±	0.022
1,2-Dimethylphenanthrene <sup>(c,d)</sup>	0.05	±	0.02
1,6-, 1,7-, 2,5-, and 2,9-Dimethylphenanthrene <sup>(c,d)</sup>	0.57	±	0.08
1,8-Dimethylphenanthrene <sup>(c,d)</sup>	0.06	±	0.02
2,6-Dimethylphenanthrene <sup>(c,d)</sup>	0.25	±	0.05
2,7-Dimethylphenanthrene <sup>(c,d)</sup>	0.23	±	0.05
3,6-Dimethylphenanthrene <sup>(c,d)</sup>	0.18	±	0.02
Fluoranthene <sup>(b,c,d,e,g,h,i)</sup>	26.9	±	1.1
Pyrene <sup>(b,c,d,e,g,h,i)</sup>	0.94	±	0.11
1-, 3-, and 7-Methylfluoranthene <sup>(c)</sup>	0.53	±	0.03
2-Methylpyrene <sup>(c,d)</sup>	0.040	±	0.008
4-Methylpyrene <sup>(c,d)</sup>	0.022	±	0.005
Benzo[ <i>c</i> ]phenanthrene <sup>(c,d,e,g,h,i,k)</sup>	1.15	±	0.21
Benz[ <i>a</i> ]anthracene <sup>(b,c,d,e,g,h,i)</sup>	0.330	±	0.037
Chrysene <sup>(b,e,k)</sup>	4.58	±	0.17
Benzo[ <i>k</i> ]fluoranthene <sup>(b,c,d,e,h,i,j,k)</sup>	0.699	±	0.058
Benzo[ <i>e</i> ]pyrene <sup>(c,d,e,g,h,i,k)</sup>	1.14	±	0.06
Benzo[ <i>a</i> ]pyrene <sup>(b,c,d,e,k)</sup>	0.053	±	0.003
Perylene <sup>(b,e)</sup>	0.059	±	0.005
Indeno[1,2,3- <i>cd</i> ]pyrene <sup>(b,c,d,e)</sup>	1.37	±	0.16
Benzo[ <i>ghi</i> ]perylene <sup>(b,c,d,e,g,h,i)</sup>	0.489	±	0.032
Dibenz[ <i>a,j</i> ]anthracene <sup>(b,c,e)</sup>	0.337	±	0.071
Benzo[ <i>b</i> ]chrysene <sup>(c,d,e)</sup>	0.094	±	0.018
Coronene <sup>(c,e)</sup>	1.32	±	0.27
PFE at 200 °C			
Naphthalene <sup>(l)</sup>	4.00	±	0.13 <sup>(f)</sup>
Biphenyl <sup>(l)</sup>	2.94	±	0.14 <sup>(f)</sup>
Fluorene <sup>(l)</sup>	2.88	±	0.19
Phenanthrene <sup>(l)</sup>	20.7	±	0.3
Anthracene <sup>(l)</sup>	0.049	±	0.002 <sup>(f)</sup>
3-Methylphenanthrene <sup>(l)</sup>	1.38	±	0.08
9-Methylphenanthrene <sup>(l)</sup>	0.587	±	0.040
Fluoranthene <sup>(l)</sup>	30.9	±	0.5
Pyrene <sup>(l)</sup>	1.52	±	0.20
Benzo[ <i>c</i> ]phenanthrene <sup>(l)</sup>	1.64	±	0.10 <sup>(f)</sup>
Benz[ <i>a</i> ]anthracene <sup>(l)</sup>	0.966	±	0.043 <sup>(f)</sup>
Chrysene <sup>(l)</sup>	5.73	±	0.11 <sup>(f)</sup>
Benzo[ <i>k</i> ]fluoranthene <sup>(l)</sup>	1.75	±	0.09 <sup>(f)</sup>
Benzo[ <i>e</i> ]pyrene <sup>(l)</sup>	2.28	±	0.18
Benzo[ <i>a</i> ]pyrene <sup>(l)</sup>	0.766	±	0.039 <sup>(f)</sup>
Perylene <sup>(l)</sup>	0.082	±	0.006
Indeno[1,2,3- <i>cd</i> ]pyrene <sup>(l)</sup>	2.13	±	0.12 <sup>(f)</sup>
Benzo[ <i>ghi</i> ]perylene <sup>(l)</sup>	1.59	±	0.09 <sup>(f)</sup>
Dibenz[ <i>a,j</i> ]anthracene <sup>(l)</sup>	0.433	±	0.008 <sup>(f)</sup>
Benzo[ <i>b</i> ]chrysene <sup>(l)</sup>	0.210	±	0.011 <sup>(f)</sup>
Coronene <sup>(l)</sup>	2.12	±	0.08 <sup>(f)</sup>

Table 3. Reference Mass Fraction Values for Selected PAHs in SRM 2975  
Based on Extraction Method and Conditions (continued)

Extraction Conditions	Mass Fractions <sup>(a)</sup> (mg/kg)
Soxhlet Extraction or PFE using temperatures between 100 °C and 200 °C	
2-Methylnaphthalene <sup>(c,d,e,l)</sup>	1.91 ± 0.32
1-Methylnaphthalene <sup>(c,d,e,l)</sup>	0.99 ± 0.11
Acenaphthene <sup>(e,l)</sup>	0.548 ± 0.029
Benzo[ghi]fluoranthene <sup>(e,l)</sup>	10.6 ± 0.3
Benzo[b]fluoranthene <sup>(d,i,k)</sup>	11.5 ± 2.8

<sup>(a)</sup> The reference values, unless otherwise footnoted, are weighted means of the mass fractions from two to six analytical methods [19]. The uncertainty listed with each value is an expanded uncertainty about the mean [19,20], with coverage factor,  $k = 2$ , calculated by combining within method variances with a between method variance [21] following the ISO/JCGM Guide [22,23].

<sup>(b)</sup> LC-FL of isomeric PAH fractions after Soxhlet extraction with DCM

<sup>(c)</sup> GC/MS (IVa) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(d)</sup> GC/MS (IVb) on 50 % phenyl-substituted methylpolysiloxane phase of same extracts as GC/MS (IVa)

<sup>(e)</sup> GC/MS (V) on non-polar extra low bleed proprietary phase after PFE with toluene at one temperature and two pressures (100 °C with 13.8 MPa and 100 °C with 20.7 MPa)

<sup>(f)</sup> The reference value is a weighted mean of average mass fractions, with one average from each of two or more analytical methods [19,20]. The expanded uncertainty is the half width of a symmetric 95 % parametric bootstrap confidence interval [24], which is consistent with the ISO/JCGM Guide [22,23]. The effective coverage factor,  $k = 2$ .

<sup>(g)</sup> GC/MS (Ia) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with toluene:methanol mixture (100 °C with 13.8 MPa)

<sup>(h)</sup> GC/MS (II) on 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(i)</sup> GC/MS (III) on 5 % phenyl-substituted methylpolysiloxane phase after Soxhlet extraction with DCM

<sup>(j)</sup> GC/MS (Ib) on 50 % phenyl-substituted methylpolysiloxane phase of selected extracts from GC/MS (Ia)

<sup>(k)</sup> GC/MS (IVc) on a smectic liquid crystalline phase of same extracts as GC/MS (IVa)

<sup>(l)</sup> GC/MS (V) on non-polar extra low bleed proprietary phase after PFE with toluene at one temperature and two pressures (200 °C with 13.8 MPa and 200 °C with 20.7 MPa)

Table 4. Reference Mass Fraction Values for Selected Nitro-PAHs in SRM 2975

Extraction Conditions	Mass Fractions <sup>(a)</sup> (mg/kg)
PFE using temperatures between 100 °C and 200 °C	
1-Nitronaphthalene <sup>(b)</sup>	0.044 ± 0.002
2-Nitronaphthalene <sup>(b)</sup>	0.118 ± 0.005
PFE at 100 °C	
9-Nitroanthracene <sup>(c)</sup>	3.07 ± 0.28
8-Nitrofluoranthene <sup>(d)</sup>	0.60 ± 0.13
4-Nitropyrene <sup>(d)</sup>	0.175 ± 0.016 <sup>(e)</sup>
6-Nitrobenzo[ <i>a</i> ]pyrene <sup>(f)</sup>	1.36 ± 0.40 <sup>(e)</sup>
1,3-Dinitropyrene <sup>(f)</sup>	1.12 ± 0.21
1,6-Dinitropyrene <sup>(f)</sup>	2.35 ± 0.51 <sup>(e)</sup>
1,8-Dinitropyrene <sup>(f)</sup>	3.10 ± 0.95
PFE at 200 °C	
9-Nitroanthracene <sup>(g)</sup>	3.94 ± 0.19 <sup>(e)</sup>

<sup>(a)</sup> The reference values, unless otherwise footnoted, are weighted means of the mass fractions from two to six analytical methods [19]. The uncertainty listed with each value is an expanded uncertainty about the mean [19,20], with coverage factor,  $k=2$ , calculated by combining within method variances with a between method variance [21] following the ISO/JCGM Guide [22,23].

<sup>(b)</sup> Two methods using GC/NCI-MS on a 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa) and four methods using GC/NCI-MS on a 50 % phenyl-substituted methylpolysiloxane phase after PFE with toluene at two temperatures and two pressures (100 °C with 13.8 MPa; 100 °C with 20.7 MPa; 200 °C with 13.8 MPa; and 200 °C with 20.7 MPa)

<sup>(c)</sup> Three methods using GC/NCI-MS on either a 5 % phenyl-substituted or 50 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa) and two methods using GC/NCI-MS on a 50 % phenyl-substituted methylpolysiloxane phase after PFE with toluene at one temperatures and two pressures (100 °C with 13.8 MPa and 100 °C with 20.7 MPa)

<sup>(d)</sup> Two methods using GC/NCI-MS on a 5 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(e)</sup> The reference value is a weighted mean of average mass fractions, with one average from each of two or more analytical methods [19,20]. The expanded uncertainty is the half width of a symmetric 95 % parametric bootstrap confidence interval [24], which is consistent with the ISO/JCGM Guide [22,23] and the effective coverage factor,  $k$ , equals 2.

<sup>(f)</sup> Three methods using GC/NCI-MS on either a 5 % phenyl-substituted or 50 % phenyl-substituted methylpolysiloxane phase after PFE with DCM (100 °C with 13.8 MPa)

<sup>(g)</sup> Two methods using GC/NCI-MS on a 50 % phenyl-substituted methylpolysiloxane phase after PFE with toluene at one temperature and two pressures (200 °C with 13.8 MPa and 200 °C with 20.7 MPa)

Table 5. Reference Mass Fraction Value for Total Extractable Mass for SRM 2975

	Mass Fraction <sup>(a,b)</sup> (%)
Total Extractable Mass	2.7 ± 0.2

<sup>(a)</sup> Extractable mass as determined from Soxhlet extraction using DCM

<sup>(b)</sup> The results are expressed as the reference value ± the expanded uncertainty. The reference value for the total extractable mass is the mean value of six measurements. The uncertainty, computed according to the CIPM approach as described in the NIST and ISO/JCGM Guide [3,22,23], is an expanded uncertainty at the 95 % level of confidence with coverage factor,  $k=2$ . The expanded uncertainty defines a range that contains the best estimate of the true value at a level of confidence of approximately 95 %.



Table 6. Reference Values for Particle-Size Characteristics for SRM 2975

Particle Measurement	Value <sup>(a)</sup>		
Mean diameter (volume distribution, MV, $\mu\text{m}$ ) <sup>(b)</sup>	31.9	±	0.6
Mean diameter (area distribution, $\mu\text{m}$ ) <sup>(c)</sup>	11.2	±	0.1
Mean diameter (number distribution, $\mu\text{m}$ ) <sup>(d)</sup>	1.62	±	0.01
Surface Area ( $\text{m}^2/\text{cm}^3$ ) <sup>(e)</sup>	0.538	±	0.006
Percentile <sup>(f)</sup>	Particle Diameter <sup>(a)</sup> ( $\mu\text{m}$ )		
95	110	±	3
90	70	±	2
80	44.9	±	0.8
70	32.4	±	0.6
60	24.8	±	0.4
50 <sup>(g)</sup>	19.4	±	0.3
40	15.2	±	0.2
30	11.7	±	0.2
20	8.5	±	0.1
10	5.3	±	0.1

<sup>(a)</sup> Each reference value is the mean value of measurements from the analysis of subsamples from four bottles. Each uncertainty, computed according to the CIPM approach as described in the ISO/JCGM Guide [3,22,23], is an expanded uncertainty at the 95 % level of confidence with coverage factor,  $k = 2$ . The expanded uncertainty defines a range that contains the best estimate of the true value at a level of confidence of approximately 95 %.

<sup>(b)</sup> The mean diameter of the volume distribution represents the center of gravity of the distribution and compensates for scattering efficiency and refractive index. This parameter is strongly influenced by coarse particles.

<sup>(c)</sup> The mean diameter of the area distribution, calculated from the volume distribution with less influence from the presence of coarse particles than the MV parameter.

<sup>(d)</sup> The mean diameter of the number distribution calculated from the volume distribution.

<sup>(e)</sup> Calculated specific surface area assuming solid, spherical particles. This is a computation and should not be interchanged with an adsorption method of surface area determination (see Table 8) as this value does not reflect porosity or topographical characteristics.

<sup>(f)</sup> The data shown is the percent of the volume that is smaller than the indicated size.

<sup>(g)</sup> Median diameter (50 % of the volume is less than 19.4  $\mu\text{m}$ )

Table 7. Information Mass Fraction Values of Selected PAHs of Relative Molecular Mass 302 in SRM 2975

	Mass Fraction <sup>(a)</sup> (mg/kg)
Dibenzo[ <i>b,e</i> ]fluoranthene	2.7
Naphtho[1,2- <i>b</i> ]fluoranthene	7.3
Naphtho[1,2- <i>k</i> ]fluoranthene and Naphtho[2,3- <i>j</i> ]fluoranthene	2.6
Naphtho[2,3- <i>b</i> ]fluoranthene	0.37
Dibenzo[ <i>a,k</i> ]fluoranthene	0.71
Dibenzo[ <i>j,l</i> ]fluoranthene	2.3
Naphtho[2,3- <i>e</i> ]pyrene	0.56
Dibenzo[ <i>e,l</i> ]pyrene	3.4

<sup>(a)</sup> These mass fraction values are provided for informational purposes only. The values have not been confirmed by an independent analytical technique as required for certification.

Table 8. Information Value for Specific Surface Area of SRM 2975

Specific Surface Area (S)<sup>(a)</sup>91 m<sup>2</sup>/g<sup>(a)</sup> Specific surface area determined by multi-point N<sub>2</sub> gas adsorption BET method.

## REFERENCES

- [1] SRM 1975; *Diesel Particulate Extract*; National Institute of Standards and Technology; U.S. Department of Commerce: Gaithersburg, MD (07 November 2000).
- [2] SRM 1650b; *Diesel Particulate Matter*; National Institute of Standards and Technology; U.S. Department of Commerce: Gaithersburg, MD (27 September 2006).
- [3] May, W.; Parris, R.; Beck II, C.; Fassett, J.; Greenberg, R.; Guenther, F.; Kramer, G.; Wise, S.; Gills, T.; Colbert, J.; Gettings, R.; MacDonald, B.; *Definition of Terms and Modes Used at NIST for Value Assignment of Reference Materials for Chemical Measurements*; NIST Special Publication 260-136 (2000); available at <http://www.nist.gov/srm/upload/SP260-136.PDF> (accessed Dec 2013).
- [4] Wright, M.E.; Klein, A.D. Jr.; Stesniak, E.J.; *A Diesel Exhaust Filter System for Industrial Diesel Forklifts*; SAE Technical Paper Series 911852: Warrendale, PA (1991).
- [5] Schantz, M.M.; Nichols, J.J.; Wise, S.A.; *Evaluation of Pressurized Fluid Extraction for the Extraction of Environmental Matrix Reference Materials*; Anal. Chem., Vol. 69, pp. 4210–4219 (1997).
- [6] Wise, S.A.; Chesler, S.N.; Hertz, H.S.; Hilpert, L.R.; May, W.E.; *Chemically-Bonded Aminosilane Stationary Phase for the High Performance Liquid Chromatographic Separation of Polynuclear Aromatic Hydrocarbons*; Anal. Chem., Vol. 49, pp. 2306–2310 (1977).
- [7] May, W.E.; Wise, S.A.; *Liquid Chromatographic Determination of Polycyclic Aromatic Hydrocarbons in Air Particulate Extracts*; Anal. Chem., Vol. 56, pp. 225–232 (1984).
- [8] Wise, S.A.; Benner, B.A. Jr.; Byrd, G.D.; Chesler, S.N.; Rebbert, R.E.; Schantz, M.M.; *Determination of Polycyclic Aromatic Hydrocarbons in a Coal Tar Standard Reference Material*; Anal. Chem., Vol. 60, pp. 887–894 (1988).
- [9] Wise, S.A.; Deissler, A.; Sander, L.C.; *Liquid Chromatographic Determination of Polycyclic Aromatic Hydrocarbon Isomers of Molecular Weight 278 and 302 in Environmental Standard Reference Materials*; Polycyclic Aromat. Compd., Vol. 3, pp. 169–184 (1993).
- [10] Scheepers, P.T.J.; Velders, D.D.; Martens, M.H.J.; Noordhoek, J.; Bos, R.P.; *Gas Chromatographic-Mass Spectrometric determination of Nitro Polycyclic Aromatic Hydrocarbons in Airborne Particulate Matter from Workplace Atmospheres Contaminated with Diesel Exhaust*; J. Chromatogr., Vol. 677, pp. 107–121 (1994).
- [11] Schantz, M.M.; McGaw, E.; Wise, S.A.; *Pressurized Liquid Extraction of Diesel and Air Particulate Standard Reference Materials: Effects of Extraction Temperature and Pressure*; Anal. Chem., Vol. 84, pp. 8222–8231 (2012).
- [12] Bergvall, C.; Westerholm, R.; *Determination of 252-302 Da and Tentative Identification of 316-376 Da Polycyclic Aromatic Hydrocarbons in Standard Reference Materials 1649a Urban Dust and 1650b and 2975 Diesel Particulate Matter by Accelerated Solvent Extraction- HPLC-GC-MS*; Anal. Bioanal. Chem., Vol. 391, pp. 2235–2248 (2008).
- [13] Masala, S.; Ahmed, T.; Bergvall, C.; Westerholm, R.; *Improved Efficiency of Extraction of Polycyclic Aromatic Hydrocarbons (PAHs) from the National Institute of Standards and Technology (NIST) Standard Reference Material Diesel Particulate Matter (SRM 2975) using Accelerated Solvent Extraction*; Anal. Bioanal. Chem., Vol. 401, pp. 3305–3315 (2011).
- [14] Cooper, L.R.; Haverland, R.L.; Vendricks, D.M.; Knisel, W.G.; *Microtrac Particle-Size Analyzer: An Alternative Particle-Size Determination Method for Sediments and Soils*; Soil Sci., Vol. 138, pp. 138–146 (1984).
- [15] Gregg, S.J.; Sing, K.S.W.; *Adsorption, Surface Area and Porosity*; 2nd ed., Academic Press, London (1982).
- [16] Brunauer, S.; Emmett, P.H.; Teller, E.; *Adsorption of Gases in Multimolecular Layers*; J. Amer. Chem. Soc., Vol. 60, pp. 309–319 (1938).
- [17] Barrett, E.P.; Joyner, L.G.; Halenda, P.P.; *The Determination of Pore Volume and Area Distributions in Porous Substances. I. Computations for Nitrogen Isotherms*; J. Amer. Chem. Soc., Vol. 73, pp. 373–380 (1951).
- [18] Hughes, T.J.; Lewtas, J.; Claxton, L.D.; *Development of a Standard Reference Material for Diesel Mutagenicity in the Salmonella Plate Incorporation Assay*; Mutat. Res., Vol. 391, pp. 243–258 (1997).
- [19] Dersimonian, R.; Laird, N.; *Meta-Analysis in Clinical Trials*; Control Clin. Trials, Vol. 7, pp. 177–188 (1986).
- [20] Rukhin, A.L.; *Weighted Means Statistics in Interlaboratory Studies*; Metrologia, Vol. 46, pp. 323–331 (2009).
- [21] Horn, R.A.; Horn, S.A.; Duncan, D.B.; *Estimating Heteroscedastic Variance in Linear Models*; J. Am. Stat. Assoc., Vol. 70, pp. 380–385 (1975).

- [22] JCGM 100:2008; *Evaluation of Measurement Data — Guide to the Expression of Uncertainty in Measurement* (GUM 1995 with Minor Corrections); Joint Committee for Guides in Metrology (2008); available at [http://www.bipm.org/utis/common/documents/jcgm/JCGM\\_100\\_2008\\_E.pdf](http://www.bipm.org/utis/common/documents/jcgm/JCGM_100_2008_E.pdf) (accessed Dec 2013); see also Taylor, B.N.; Kuyatt, C.E.; *Guidelines for Evaluating and Expressing the Uncertainty of NIST Measurement Results*; NIST Technical Note 1297; U.S. Government Printing Office: Washington, DC (1994); available at <http://www.nist.gov/phylab/pubs/index.cfm> (accessed Dec 2013).
- [23] JCGM 101:2008; *Evaluation of Measurement Data – Supplement 1 to the “Guide to the Expression of Uncertainty in Measurement” – Propagation of Distributions using a Monte Carlo method*; JCGM (2008); available at [http://www.bipm.org/utis/common/documents/jcgm/JCGM\\_101\\_2008\\_E.pdf](http://www.bipm.org/utis/common/documents/jcgm/JCGM_101_2008_E.pdf) (accessed Dec 2013).
- [24] Efron, B.; Tibshirani, R.J.; *An Introduction to the Bootstrap*; Chapman & Hall (1993).

**Certificate Revision History:** 03 December 2013 (Update certified and reference values for PAHs and nitro-PAHs to address the effect of temperature used for extraction; extension of certification period; editorial changes); 19 March 2009 (Addition of reference and information values and update of expiration date); 10 April 2008 (Update of expiration date and editorial changes); 07 November 2000 (Original certificate date).

*Users of this SRM should ensure that the Certificate of Analysis in their possession is current. This can be accomplished by contacting the SRM Program: telephone (301) 975-2200; fax (301) 948-3730; e-mail [srminfo@nist.gov](mailto:srminfo@nist.gov); or via the Internet at <http://www.nist.gov/srm>.*