

JOINT RESEARCH CENTRE  
Directorate F – Health, Consumers and Reference Materials

# REFERENCE MATERIAL CERTIFICATE

## ERM® - FD305

SILICA PARTICLES IN AQUEOUS SOLUTION		
	Certified value <sup>3)</sup>	Uncertainty <sup>4)</sup>
Mean electrophoretic mobility <sup>1)</sup> ( $10^{-8} \text{ m}^2 \text{ V}^{-1} \cdot \text{s}^{-1}$ )	-4.5	0.4
Mean zeta potential (mV) <sup>2)</sup>	-58	5
<p><sup>1)</sup> As obtained with electrophoretic light scattering (ELS) at a sample temperature of 25 °C and by applying ISO 13099-2:2012.</p> <p><sup>2)</sup> As calculated from the certified electrophoretic mobility value using the Smoluchowski approximation <math>f(\kappa a) = 1.5</math> (for a temperature of 25 °C, a dynamic viscosity of 0.89 mPa·s and a value of 78.4 for the dielectric constant of water).</p> <p><sup>3)</sup> Certified values are values that fulfil the highest standards of accuracy and represent the unweighted mean value of the means of accepted sets of data; each set being obtained in a different laboratory or with a different method of determination. The certified value and its uncertainty are traceable to the International System of units (SI).</p> <p><sup>4)</sup> The uncertainty of the certified value is the expanded uncertainty with a coverage factor <math>k = 2</math> corresponding to a level of confidence of about 95 % estimated in accordance with ISO 17034:2016 and ISO Guide 35:2017.</p>		

This certificate is valid for one year after purchase.

Sales date:

The minimum amount of sample to be used is 400 µL.

**Note:** This material is a joint production of the European Commission's Joint Research Centre (JRC) and the US National Institute of Standards and Technology (NIST). ERM-FD305 is identical to NIST SRM 1992.

Geel, October 2020

Signed: 

Dr Robert Koeber  
Head of Unit Reference Materials  
European Commission, Joint Research Centre  
Directorate F – Health, Consumers and Reference Materials  
Retieseweg 111  
B-2440 Geel, Belgium

Additional Material Information	
	Value (at 25°C)
Conductivity <sup>1)</sup>	0.39 mS/cm
1) Mean of the conductivity measurement results reported by the laboratories.	

**Note :** Zeta potential values as obtained by particle tracking analysis (PTA) and electrical asymmetrical flow-field flow fractionation (EAF4) are available in the certification report of the material.

## DESCRIPTION OF THE MATERIAL

ERM-FD305 consists of a suspension of silica particles at a mass fraction of 1.5 g/kg in a borate buffer at pH 8.9. The material is available in 5 mL pre-scored ampoules containing approximately 5 mL of suspension.

## ANALYTICAL METHODS USED FOR CHARACTERISATION

Electrophoretic light scattering (ELS)

Particle tracking analysis (PTA)

Electrical asymmetrical flow-field flow fractionation (EA4F)

## PARTICIPANTS

Anton Paar GmbH, Graz, AT

European Commission, Joint Research Centre, Directorate F – Health, Consumers and Reference Materials, Geel, BE

(measurements under the scope of ISO/IEC 17025 accreditation BELAC 268-TEST)

Fraunhofer Institut für Keramische Technologie and Systeme (IKTS), Dresden, DE

(Measurements under the scope of ISO/IEC 17025 accreditation DAkkS D-PL-11140-15-00)

Horiba, Palaiseau, FR

LGC Ltd, Teddington, UK

Malvern Panalytical Ltd, Malvern, UK

Malvern Panalytical Inc., Westborough, USA

Microtrac Inc., Montgomeryville, USA

Moscow Institute of Physics and Technology, Moscow, RU

National Institute of Standards and Technology (NIST), Gaithersburg, USA

National Measurement Institute Australia (NMIA), West Lindfield, AU

National Metrology Institute of Japan (NMIJ), Tsukuba, JPN

National Physical Laboratory (NPL), Teddington, UK

Otsuka Electronics, Osaka, JPN

Postnova Analytics GmbH, Landsberg, DE

Wyatt Technology, Santa Barbara, USA

## SAFETY INFORMATION

The material should be handled with care. The material contains a fraction of nanoparticles (size range of 1 nm to 100 nm) that could have an impact on environment and human health. Any spillage of the suspension should be handled according to the standard laboratory safety precautions.

For further details refer to the safety data sheet.

## INTENDED USE

The main purpose of this material is to assess performance of instruments and/or methods that are used for measuring zeta potential and electrophoretic mobility. As with any reference material, it can be used for establishing control charts or validation studies.

## INSTRUCTIONS FOR USE

Before opening, the ampoule should be gently inverted several times to ensure the homogeneity of the suspension and to re-suspend any settled particles. Remove any suspension that remains in the upper part (conical top) of the ampoule by gently flicking the conical part with the forefinger while tilting the ampoule. The ampoule is pre-scored and can be opened by applying gentle pressure with one's thumb to snap off the conical part. The content of the ampoule should be used the same day as opened and should be gently homogenised before every measurement without introducing air bubbles.

ELS method: Aliquot of ERM-FD305 should be measured as received, i.e. without dilution. The measurement temperature shall be  $25\text{ °C} \pm 0.2\text{ °C}$ . Values to be used for the viscosity and refractive index of the dispersing medium (water) at  $25\text{ °C}$  are  $0.89\text{ mPa}\cdot\text{s}$  and  $1.332$ , respectively. The value of the viscosity must be adjusted when tests are not performed at  $25\text{ °C}$ .

For general information on handling of reference materials, please see ERM Application Note 6, available on <https://crm.jrc.ec.europa.eu/e/132/User-support-Application-Notes>.

## STORAGE

The materials should be stored at ambient temperature ( $20 \pm 5$ ) °C. Ampoules must not be allowed to freeze, as this will irreversibly compromise the integrity of the material.

For more information regarding the shelf life of reference materials please see ERM Application Note 7, available on <https://crm.jrc.ec.europa.eu/e/132/User-support-Application-Notes>.

Please note that the stability of opened samples has not been tested and repeated use of the material occurs under the responsibility of the user. The European Commission cannot be held responsible for changes that happen during storage of the material at the customer's premises, especially of opened samples.

## LEGAL NOTICE

Neither the European Commission, its contractors nor any person acting on their behalf:

(a) make any warranty or representation, express or implied, that the use of any information, material, apparatus, method or process disclosed in this document does not infringe any privately owned intellectual property rights; or

(b) assume any liability with respect to, or for damages resulting from, the use of any information, material, apparatus, method or process disclosed in this document save for loss or damage arising solely and directly from the negligence of the Joint Research Centre of the European Commission.

## NOTE

A detailed certification report is available at <https://crm.jrc.ec.europa.eu/>.

A paper copy can be obtained from the Joint Research Centre, Directorate F – Health, Consumers and Reference Materials on request.



---

European Commission – Joint Research Centre  
Directorate F – Health, Consumers and Reference Materials  
Retieseweg 111, B - 2440 Geel (Belgium)  
Telephone: +32-(0)14-571.705 - Fax: +32-(0)14-590.406  
[jrc-rm-distribution@ec.europa.eu](mailto:jrc-rm-distribution@ec.europa.eu)