

Congener Calibration Set

Art. ID ACCU-C-CSA-SET
Unit 5 x 1 mL
Deliverydetails De Minimis

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-C-CS-01	2,2',3,3',4,4',5,5',6-N onachlorobiphenyl (PCB 206)	[40186-72-9]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4,4',5,5'-Oct achlorobiphenyl (PCB 19 4)	[35694-08-7]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4,4',5,6-Octa chlorobiphenyl (PCB 195)	[52663-78-2]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4',5,5',6-Oct achlorobiphenyl (IUPAC PCB 199, PCB 201)	[52663-75-9]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4,5,6-Heptach lorobiphenyl (PCB 173)	[68194-16-1]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4,5,6'-Heptac hlorobiphenyl (PCB 174)	[38411-25-5]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4',5,6-Heptac hlorobiphenyl (PCB 177)	[52663-70-4]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',4-Pentachloro biphenyl (PCB 82)	[52663-62-4]	10	µg/mL		
ACCU-C-CS-01	2,2',3,3',5,6,6'-Heptac hlorobiphenyl (PCB 179)	[52663-64-6]	10	µg/mL		
ACCU-C-CS-01	2,2',3,4,4',5,5',6-Octa chlorobiphenyl (PCB 203)	[52663-76-0]	10	µg/mL		
ACCU-C-CS-01	2,2',3,4,4',5,5'-Heptac hlorobiphenyl (PCB 180)	[35065-29-3]	10	µg/mL		
ACCU-C-CS-01	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	10	µg/mL		
ACCU-C-CS-01	2,2',3,4',5,5',6-Heptac hlorobiphenyl (PCB 187)	[52663-68-0]	10	µg/mL		

ACCU-C-CS-01	2,2',3,4',5,5'-Hexachlorobiphenyl (PCB 146)	[51908-16-8]	10	µg/mL
ACCU-C-CS-01	2,2',3,4',5,6-Hexachlorobiphenyl (PCB 147)	[68194-13-8]	10	µg/mL
ACCU-C-CS-01	2,2',3,4,5'-Pentachlorobiphenyl (PCB 87)	[38380-02-8]	10	µg/mL
ACCU-C-CS-01	2,2',3,5'-Tetrachlorobiphenyl (PCB 44)	[41464-39-5]	10	µg/mL
ACCU-C-CS-01	2,2',3-Trichlorobiphenyl I (PCB 16)	[38444-78-9]	10	µg/mL
ACCU-C-CS-01	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	10	µg/mL
ACCU-C-CS-01	2,2',4,4',5-Pentachlorobiphenyl (PCB 99)	[38380-01-7]	10	µg/mL
ACCU-C-CS-01	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	10	µg/mL
ACCU-C-CS-01	2,2',5-Trichlorobiphenyl I (PCB 18)	[37680-65-2]	10	µg/mL
ACCU-C-CS-01	2,2',6-Trichlorobiphenyl I (PCB 19)	[38444-73-4]	10	µg/mL
ACCU-C-CS-01	2,2'-Dichlorobiphenyl (PCB 4)	[13029-08-8]	10	µg/mL
ACCU-C-CS-01	2,3,3',4',6-Pentachlorobiphenyl (PCB 110)	[38380-03-9]	10	µg/mL
ACCU-C-CS-01	2,3,3',4'-Tetrachlorobiphenyl (PCB 56)	[41464-43-1]	10	µg/mL
ACCU-C-CS-01	2,3',4,4'-Tetrachlorobiphenyl (PCB 66)	[32598-10-0]	10	µg/mL
ACCU-C-CS-01	2,3',4,5-Tetrachlorobiphenyl (PCB 67)	[73575-53-8]	10	µg/mL
ACCU-C-CS-01	2,3',4',6-Tetrachlorobiphenyl (PCB 71)	[41464-46-4]	10	µg/mL
ACCU-C-CS-01	2,3,4'-Trichlorobiphenyl I (PCB 22)	[38444-85-8]	10	µg/mL
ACCU-C-CS-01	2,3',4-Trichlorobiphenyl I (PCB 25)	[55712-37-3]	10	µg/mL
ACCU-C-CS-01	2,3'-Dichlorobiphenyl (PCB 6)	[25569-80-6]	10	µg/mL

ACCU-C-CS-01	2,4,4',5-Tetrachlorobiphenyl (PCB 74)	[32690-93-0]	10	µg/mL
ACCU-C-CS-01	2,4,4'-Trichlorobiphenyl (PCB 28)	[7012-37-5]	10	µg/mL
ACCU-C-CS-01	2,4'-Dichlorobiphenyl (PCB 8)	[34883-43-7]	10	µg/mL
ACCU-C-CS-01	2,5-Dichlorobiphenyl (PCB 9)	[34883-39-1]	10	µg/mL
ACCU-C-CS-01	2-Chlorobiphenyl (PCB 1)	[2051-60-7]	10	µg/mL
ACCU-C-CS-01	3-Chlorobiphenyl (PCB 2)	[2051-61-8]	10	µg/mL
ACCU-C-CS-01	4-Chlorobiphenyl (PCB 3)	[2051-62-9]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,4',5,6-Octachlorobiphenyl (PCB 196)	[42740-50-1]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,4',5-Heptachlorobiphenyl (PCB 170)	[35065-30-6]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,4',6-Octachlorobiphenyl (PCB 197)	[33091-17-7]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,4',6-Heptachlorobiphenyl (PCB 171)	[52663-71-5]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,5,5'-Heptachlorobiphenyl (PCB 172)	[52663-74-8]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,6-Hexachlorobiphenyl (PCB 131)	[61798-70-7]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',4,6'-Hexachlorobiphenyl (PCB 132)	[38380-05-1]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',5,5',6-Heptachlorobiphenyl (PCB 178)	[52663-67-9]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',5,6'-Hexachlorobiphenyl (PCB 135)	[52744-13-5]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',5-Pentachlorobiphenyl (PCB 83)	[60145-20-2]	10	µg/mL
ACCU-C-CS-02	2,2',3,3',6-Pentachlorobiphenyl (PCB 84)	[52663-60-2]	10	µg/mL

ACCU-C-CS-02	2,2',3,4,4',5',6-Heptachlorobiphenyl (PCB 183)	[52663-69-1]	10	µg/mL
ACCU-C-CS-02	2,2',3,4,5,5'-Hexachlorobiphenyl (PCB 141)	[52712-04-6]	10	µg/mL
ACCU-C-CS-02	2,2',3,4',5',6-Hexachlorobiphenyl (PCB 149)	[38380-04-0]	10	µg/mL
ACCU-C-CS-02	2,2',3,4-Tetrachlorobiphenyl (PCB 41)	[52663-59-9]	10	µg/mL
ACCU-C-CS-02	2,2',3,6-Tetrachlorobiphenyl (PCB 45)	[70362-45-7]	10	µg/mL
ACCU-C-CS-02	2,2',3,6'-Tetrachlorobiphenyl (PCB 46)	[41464-47-5]	10	µg/mL
ACCU-C-CS-02	2,2',4,5',6-Pentachlorobiphenyl (PCB 103)	[60145-21-3]	10	µg/mL
ACCU-C-CS-02	2,2',4,5-Tetrachlorobiphenyl (PCB 48)	[70362-47-9]	10	µg/mL
ACCU-C-CS-02	2,2',4-Trichlorobiphenyl (PCB 17)	[37680-66-3]	10	µg/mL
ACCU-C-CS-02	2,3,3',4,4',5,5',6-Octachlorobiphenyl (PCB 205)	[74472-53-0]	10	µg/mL
ACCU-C-CS-02	2,3,3',4',5,5',6-Heptachlorobiphenyl (PCB 193)	[69782-91-8]	10	µg/mL
ACCU-C-CS-02	2,3,3',4',5',6-Hexachlorobiphenyl (PCB 164)	[74472-45-0]	10	µg/mL
ACCU-C-CS-02	2,3,3',4',5-Pentachlorobiphenyl (PCB 107)	[70424-68-9]	10	µg/mL
ACCU-C-CS-02	2,3,4,4',6-Pentachlorobiphenyl (PCB 115)	[74472-38-1]	10	µg/mL
ACCU-C-CS-02	2,3,4,4'-Tetrachlorobiphenyl (PCB 60)	[33025-41-1]	10	µg/mL
ACCU-C-CS-02	2,3',4',5-Tetrachlorobiphenyl (PCB 70)	[32598-11-1]	10	µg/mL
ACCU-C-CS-02	2,3',5-Trichlorobiphenyl (PCB 26)	[38444-81-4]	10	µg/mL
ACCU-C-CS-02	2,3,6-Trichlorobiphenyl (PCB 24)	[55702-45-9]	10	µg/mL
ACCU-C-CS-02	2,3-Dichlorobiphenyl (P	[16605-91-7]	10	µg/mL

	CB 5)			
ACCU-C-CS-02	2,4',5-Trichlorobipheny I (PCB 31)	[16606-02-3]	10	µg/mL
ACCU-C-CS-02	2,4',6-Trichlorobipheny I (PCB 32)	[38444-77-8]	10	µg/mL
ACCU-C-CS-02	2,4-Dichlorobiphenyl (P CB 7)	[33284-50-3]	10	µg/mL
ACCU-C-CS-02	2,6-Dichlorobiphenyl (P CB 10)	[33146-45-1]	10	µg/mL
ACCU-C-CS-02	3,4,4'-Trichlorobipheny I (PCB 37)	[38444-90-5]	10	µg/mL
ACCU-C-CS-02	PCB No. 95	[38379-99-6]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',4,4',5,5',6,6 '-Decachlorobiphenyl (P CB 209)	[2051-24-3]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',4,4',5,6,6'-N onachlorobiphenyl (PCB 207)	[52663-79-3]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',4,4'-Hexachlo robiphenyl (PCB 128)	[38380-07-3]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',4,5,5',6,6'-N onachlorobiphenyl (PCB 208)	[52663-77-1]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',5,6-Hexachlor obiphenyl (PCB 134)	[52704-70-8]	10	µg/mL
ACCU-C-CS-03	2,2',3,3',6,6'-Hexachlo robiphenyl (PCB 136)	[38411-22-2]	10	µg/mL
ACCU-C-CS-03	2,2',3,3'-Tetrachlorobi phenyl (PCB 40)	[38444-93-8]	10	µg/mL
ACCU-C-CS-03	2,2',3,4,5',6-Hexachlor obiphenyl (PCB 144)	[68194-14-9]	10	µg/mL
ACCU-C-CS-03	2,2',3,4'-Tetrachlorobi phenyl (PCB 42)	[36559-22-5]	10	µg/mL
ACCU-C-CS-03	2,2',3,5,5',6-Hexachlor obiphenyl (PCB 151)	[52663-63-5]	10	µg/mL
ACCU-C-CS-03	2,2',3,5,5'-Pentachloro biphenyl (PCB 92)	[52663-61-3]	10	µg/mL
ACCU-C-CS-03	2,2',3,5,6-Pentachlorob	[73575-56-1]	10	µg/mL

	iphenyl (PCB 93)			
ACCU-C-CS-03	2,2',4,4'-Tetrachlorobi phenyl (PCB 47)	[2437-79-8]	10	µg/mL
ACCU-C-CS-03	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	10	µg/mL
ACCU-C-CS-03	2,3,3',4,4',5,6-Heptach lorobiphenyl (PCB 190)	[41411-64-7]	10	µg/mL
ACCU-C-CS-03	2,3,3',4,4',5',6-Heptac hlorobiphenyl (PCB 191)	[74472-50-7]	10	µg/mL
ACCU-C-CS-03	2,3,3',4,4',5'-Hexachlo robiphenyl (PCB 157)	[69782-90-7]	10	µg/mL
ACCU-C-CS-03	2,3,3',4,4',6-Hexachlor obiphenyl (PCB 158)	[74472-42-7]	10	µg/mL
ACCU-C-CS-03	2,3,3',4,4'-Pentachloro biphenyl (PCB 105)	[32598-14-4]	10	µg/mL
ACCU-C-CS-03	2,3,3'-Trichlorobipheny l (PCB 20)	[38444-84-7]	10	µg/mL
ACCU-C-CS-03	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	10	µg/mL
ACCU-C-CS-03	2,3',4,4',6-Pentachloro biphenyl (PCB 119)	[56558-17-9]	10	µg/mL
ACCU-C-CS-03	2,3',4,6-Tetrachlorobip henyl (PCB 69)	[60233-24-1]	10	µg/mL
ACCU-C-CS-03	2',3,5-Trichlorobipheny l (PCB 34)	[37680-68-5]	10	µg/mL
ACCU-C-CS-03	2,3',6-Trichlorobipheny l (PCB 27)	[38444-76-7]	10	µg/mL
ACCU-C-CS-03	2,4,5-Trichlorobiphenyl (PCB 29)	[15862-07-4]	10	µg/mL
ACCU-C-CS-03	4,4'-Dichlorobiphenyl (PCB 15)	[2050-68-2]	10	µg/mL
ACCU-C-CS-04	2,2',3,3',4,5,6,6'-Octa chlorobiphenyl (IUPAC P CB 200, PCB 199)	[52663-73-7]	10	µg/mL
ACCU-C-CS-04	2,2',3,3',4,5',6,6'-Oct achlorobiphenyl (IUPAC PCB 201, PCB 200)	[40186-71-8]	10	µg/mL
ACCU-C-CS-04	2,2',3,3',4,5',6-Heptac	[40186-70-7]	10	µg/mL

	hlorobiphenyl (PCB 175)			
ACCU-C-CS-04	2,2',3,3',4,5'-Hexachlorobiphenyl (PCB 130)	[52663-66-8]	10	µg/mL
ACCU-C-CS-04	2,2',3,3',5,5',6,6'-Octachlorobiphenyl (PCB 202)	[2136-99-4]	10	µg/mL
ACCU-C-CS-04	2,2',3,4',5-Pentachlorobiphenyl (PCB 90)	[68194-07-0]	10	µg/mL
ACCU-C-CS-04	2,2',4,4',5,6'-Hexachlorobiphenyl (PCB 154)	[60145-22-4]	10	µg/mL
ACCU-C-CS-04	2,2',4,4',6-Pentachlorobiphenyl (PCB 100)	[39485-83-1]	10	µg/mL
ACCU-C-CS-04	2,2',4,6'-Tetrachlorobiphenyl (PCB 51)	[68194-04-7]	10	µg/mL
ACCU-C-CS-04	2,2',5,6'-Tetrachlorobiphenyl (PCB 53)	[41464-41-9]	10	µg/mL
ACCU-C-CS-04	2,2',6,6'-Tetrachlorobiphenyl (PCB 54)	[15968-05-5]	10	µg/mL
ACCU-C-CS-04	2,3,3',4',5,6-Hexachlorobiphenyl (PCB 163)	[74472-44-9]	10	µg/mL
ACCU-C-CS-04	2',3,3',4,5-Pentachlorobiphenyl (PCB 122)	[76842-07-4]	10	µg/mL
ACCU-C-CS-04	2,3,3',5,5',6-Hexachlorobiphenyl (PCB 165)	[74472-46-1]	10	µg/mL
ACCU-C-CS-04	2',3,4,5,5'-Pentachlorobiphenyl (PCB 124)	[70424-70-3]	10	µg/mL
ACCU-C-CS-04	2,3,4',5,6-Pentachlorobiphenyl (PCB 117)	[68194-11-6]	10	µg/mL
ACCU-C-CS-04	2,3',5',6-Tetrachlorobiphenyl (PCB 73)	[74338-23-1]	10	µg/mL
ACCU-C-CS-04	2,4,4',6-Tetrachlorobiphenyl (PCB 75)	[32598-12-2]	10	µg/mL
ACCU-C-CS-04	3,3',4-Trichlorobiphenyl (PCB 35)	[37680-69-6]	10	µg/mL
ACCU-C-CS-04	3,4,4',5-Tetrachlorobiphenyl (PCB 81)	[70362-50-4]	10	µg/mL
ACCU-C-CS-04	3,4'-Dichlorobiphenyl (PCB 13)	[2974-90-5]	10	µg/mL

ACCU-C-CS-04	3,5-Dichlorobiphenyl (P CB 14)	[34883-41-5]	10	µg/mL
ACCU-C-CS-05	2,2',3,3',4,5-Hexachlor obiphenyl (PCB 129)	[55215-18-4]	10	µg/mL
ACCU-C-CS-05	2,2',3,3',4,6,6'-Heptac hlorobiphenyl (PCB 176)	[52663-65-7]	10	µg/mL
ACCU-C-CS-05	2,2',3,4,4',5-Hexachlor obiphenyl (PCB 136)	[35694-06-5]	10	µg/mL
ACCU-C-CS-05	2,2',3,4,4'-Pentachloro biphenyl (PCB 85)	[65510-45-4]	10	µg/mL
ACCU-C-CS-05	2,2',3,4,5,5',6-Heptach lorobiphenyl (PCB 185)	[52712-05-7]	10	µg/mL
ACCU-C-CS-05	2,2',3,4',5'-Pentachlor obiphenyl	[41464-51-1]	10	µg/mL
ACCU-C-CS-05	2,2',3,4',6-Pentachloro biphenyl (PCB 91)	[68194-05-8]	10	µg/mL
ACCU-C-CS-05	2,2',4,5'-Tetrachlorobi phenyl (PCB 49)	[41464-40-8]	10	µg/mL
ACCU-C-CS-05	2,2',4,6,6'-Pentachloro biphenyl (PCB 104)	[56558-16-8]	10	µg/mL
ACCU-C-CS-05	2,3,3',4,4',5,5'-Heptac hlorobiphenyl (PCB 189)	[39635-31-9]	10	µg/mL
ACCU-C-CS-05	2,3,3',4,4',5-Hexachlor obiphenyl (PCB 156)	[38380-08-4]	10	µg/mL
ACCU-C-CS-05	2,3,3',6-Tetrachlorobip henyl (PCB 59)	[74472-33-6]	10	µg/mL
ACCU-C-CS-05	2,3',4,4',5,5'-Hexachlo robiphenyl (PCB 167)	[52663-72-6]	10	µg/mL
ACCU-C-CS-05	2,3,4,4',5-Pentachlorob iphenyl (PCB 114)	[74472-37-0]	10	µg/mL
ACCU-C-CS-05	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	10	µg/mL
ACCU-C-CS-05	2,3,4',5-Tetrachlorobip henyl (PCB 63)	[74472-34-7]	10	µg/mL
ACCU-C-CS-05	2,3,4',6-Tetrachlorobip henyl (PCB 64)	[52663-58-8]	10	µg/mL
ACCU-C-CS-05	2',3,4-Trichlorobipheny l (PCB 33)	[38444-86-9]	10	µg/mL

ACCU-C-CS-05	3,3',4,4'-Tetrachlorobi phenyl (PCB 77)	[32598-13-3]	10	µg/mL
ACCU-C-CS-05	3,4-Dichlorobiphenyl (P CB 12)	[2974-92-7]	10	µg/mL