

Labmix24 GmbH Kesseldorfer Rott 24 46499 Hamminkeln Germany
 Tel:
 +49

 Fax:
 +49

 Web:
 www

 E-Mail:
 info

+49 (0) 2852 96064 00 +49 (0) 2852 96064 24 www.labmix24.com info@labmix24.com

Water Conductivity: USP General Chapter 645 (On-Demand)

Art. ID

Unit

USP-CM-645-01 Each

Deliverydetails No Dangerous Good

Description

USP General Chapter <645> Water Conductivity has been official since USP 23 (1996), and it has been adopted by other pharmacopoeias as the primary method for the determination of inorganic/ionic impurities in Purified Water and WFI. USP <645> Water Conductivity is one of the two primary chemical limit tests (USP <643> Total Organic Carbon is the other) to determine that there is sufficient control/reduction of chemicals in the water purification system. In addition, electrical conductivity is also used as an analytical on-line, real-time process control tool for the entire water purification process and other manufacturing processes. This webinar will explain the purpose of this test, instrumentation and calibration requirements, on-line and off-line measurements, temperature compensation, and test limit requirements for various pharmaceutical waters.