

XRF - Setting Up Sample

Art. ID XRF-AS1
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium oxide	[1344-28-1]	15,8	%		
	Arsenic trioxide (As ₂ O ₃)		0,05	%		
	Boron trioxide		3,22	%		
	Cadmium oxide (CdO)		0,39	%		
	Calcium Oxide (CaO)	[1305-78-8]	0,83	%		
	Cesium oxide (CsO)		0,04	%		
	Chromium(III) oxide	[1308-38-9]	0,15	%		
	Fluorine		0,17	%		
	Germanium dioxide (GeO ₂)	[1310-53-8]	0,08	%		
	Indium(III) oxide	[1312-43-2]	0,04	%		
	Iron (III) oxide		1,16	%		
	Lead(II) oxide	[1317-36-8]	0,5	%		
	Magnesium oxide (MgO)	[1309-48-4]	3,2	%		
	Magnesium oxide (MgO)	[1309-48-4]	20,3	%		
	Phosphorus Pentoxide (P ₂ O ₅)	[1314-56-3]	0,58	%		
	Potassium Oxide		2,16	%		
	Rubidium oxide		0,04	%		
	Silicon dioxide (SiO ₂)	[7631-86-9]	38,9	%		
	Sodium oxide	[1313-59-3]	0,13	%		
	Strontium oxide	[1314-11-0]	0,71	%		
	Tellurium dioxide	[7446-07-3]	0,04	%		
	Thorium dioxide		0,04	%		
	Titanium dioxide (TiO ₂)	[13463-67-7]	3,9	%		
	Triuranium octoxide		0,01	%		
	Vanadium oxide (V ₂ O ₅)		0,01	%		
	Zinc oxide (ZnO)	[1314-13-2]	7,4	%		
	Zirconium dioxide (ZrO ₂)		0,15	%		

