

Spiking Std Set Method 200,7 (1991)

Art. ID ACCU-M-200.7-SP-R-SET
Unit 5 x 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-M-200.7-SP-01-R	Boron (B)	[7440-42-8]	400	µg/mL		
ACCU-M-200.7-SP-01-R	Molybdenum (Mo)	[7439-98-7]	200	µg/mL		
ACCU-M-200.7-SP-01-R	Phosphorus (P)	[7723-14-0]	400	µg/mL		
ACCU-M-200.7-SP-01-R	Silicon (Si)	[7440-21-3]	2000	µg/mL		
ACCU-M-200.7-SP-02-R	Calcium (Ca)	[7440-70-2]	10000	µg/mL		
ACCU-M-200.7-SP-02-R	Magnesium (Mg)	[7439-95-4]	10000	µg/mL		
ACCU-M-200.7-SP-02-R	Potassium (K)	[7440-09-7]	10000	µg/mL		
ACCU-M-200.7-SP-02-R	Sodium (Na)	[7440-23-5]	10000	µg/mL		
ACCU-M-200.7-SP-03	Aluminium (Al)	[7429-90-5]	2000	µg/mL		
ACCU-M-200.7-SP-03	Barium (Ba)	[7440-39-3]	2000	µg/mL		
ACCU-M-200.7-SP-03	Beryllium (Be)	[7440-41-7]	50	µg/mL		
ACCU-M-200.7-SP-03	Chromium (Cr)	[7440-47-3]	200	µg/mL		
ACCU-M-200.7-SP-03	Cobalt (Co)	[7440-48-4]	500	µg/mL		
ACCU-M-200.7-SP-03	Copper (Cu)	[7440-50-8]	250	µg/mL		
ACCU-M-200.7-SP-03	Iron (Fe)	[7439-89-6]	1000	µg/mL		
ACCU-M-200.7-SP-03	Manganese (Mn)	[7439-96-5]	500	µg/mL		
ACCU-M-200.7-SP-03	Nickel (Ni)	[7440-02-0]	500	µg/mL		
ACCU-M-200.7-SP-03	Silver (Ag)	[7440-22-4]	50	µg/mL		
ACCU-M-200.7-SP-03	Vanadium (V)	[7440-62-2]	500	µg/mL		
ACCU-M-200.7-SP-03	Zinc (Zn)	[7440-66-6]	500	µg/mL		
ACCU-M-200.7-SP-04-R	Antimony (Sb)	[7440-36-0]	200	µg/mL		
ACCU-M-200.7-SP-05-R	Arsenic (As)	[7440-38-2]	200	µg/mL		
ACCU-M-200.7-SP-05-R	Cadmium (Cd)	[7440-43-9]	100	µg/mL		
ACCU-M-200.7-SP-05-R	Lead (Pb)	[7439-92-1]	200	µg/mL		
ACCU-M-200.7-SP-05-R	Selenium (Se)	[7782-49-2]	400	µg/mL		
ACCU-M-200.7-SP-05-R	Thallium (Tl)	[7440-28-0]	400	µg/mL		