

Apo-Direct TUNEL Assay Kit The APO-DIRECT Kit is a single-step staining method for labeling DNA breaks to detect apoptotic cells by flow cytometry.

Art. ID	SAF-APT110
Unit	EA
Deliverydetails	No Dangerous Good

Description

The CHEMICON APO-DIRECT Kit is a single-step staining method for labeling DNA breaks to detect apoptotic cells by flow cytometry (Li et al. 1995). The kit contains all the reagents required for measuring apoptosis in cells including positive and negative control cells, washing, reaction, and rinsing buffers, and propidium iodide/RNase A solution for counter staining the total DNA.