

## **Seronorm™ Urine L-2**

Art. ID	SERO-208605
Unit	10 x 8 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Urine control from human urine with stable creatinine /// Product documentation with assayed values for the most frequently used instruments on the market /// Two clinically relevant levels - each level available separately /// Liquid - ready to use from the fridge /// 31 days stability in opened vials /// Please refer to the package insert of currently available lots for specific analyte and stability claims

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Albumin/Creatinin quota					
	Amylase					
	Calcium (Ca)	[7440-70-2]				
	Chloride (Cl)					
	Cortisol, free					
	Creatinine	[60-27-5]				
	Glucose					
	Magnesium (Mg)	[7439-95-4]				
	Microalbumin					
	Osmolality					
	pH					
	Phosphorus (P)	[7723-14-0]				
	Potassium (K)	[7440-09-7]				
	Pregnancy					
	Protein, total					
	Sodium (Na)	[7440-23-5]				
	Specific Gravity					
	Urea	[57-13-6]				
	Uric Acid	[69-93-2]				