

Aflatoxins B/G, Fumonisin, Ochratoxin A and Zearalenone in maize flour

Art. ID AOK-PT-MA5600/YFOZA
Unit 80 g
Deliverydetails /not restricted

Description

Order deadline: 02.05.2025 /// Each test material may contain one or more substances from the table below. The proposed concentrations are indicative, especially for incurred matrix (naturally contaminated) with multiple contamination. In case 2 test materials are provided, one could be blank. It is not requested to research all the molecules of each group.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aflatoxin B1	[1162-65-8]	<50	µg/kg		
	Aflatoxin B2	[7220-81-7]	<50	µg/kg		
	Aflatoxin G1	[1165-39-5]	<50	µg/kg		
	Aflatoxin G2	[7241-98-7]	<50	µg/kg		
	Fumonisin B1	[116355-83-0]	<6	µg/kg		
	Fumonisin B2	[116355-84-1]	<6	µg/kg		
	Ochratoxin A	[303-47-9]	<30	µg/kg		
	Zearalenone	[17924-92-4]	<500	µg/kg		