

Iron base, Steel - Ti stabilised Stainless, CRM, Chips 100 g

Art. ID BAS-BCS-CRM 465/2
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0196	%		
	Carbon (C)	[7440-44-0]	0,04	%		
	Cobalt (Co)	[7440-48-4]	0,1711	%		
	Chromium (Cr)	[7440-47-3]	17,19	%		
	Copper (Cu)	[7440-50-8]	0,553	%		
	Manganese (Mn)	[7439-96-5]	1,613	%		
	Molybdenum (Mo)	[7439-98-7]	0,456	%		
	Nickel (Ni)	[7440-02-0]	9,15	%		
	Nickel (Ni)	[7440-02-0]	0,0141	%		
	Phosphorus (P)	[7723-14-0]	0,0349	%		
	Silicon (Si)	[7440-21-3]	0,407	%		
	Silicon (Si)	[7440-21-3]	0,0212	%		
	Tin (Sn)	[7440-31-5]	0,0113	%		
	Titanium (Ti)	[7440-32-6]	0,319	%		
	Vanadium (V)	[7440-62-2]	0,0857	%		