

QuiKResponse - TCLP Metals in Soil, priority pollutant (PriorityPollutnT TM)

| | |
|-----------------|-----------------------------------|
| Art. ID | ERA-544QR |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 105 g soil standard in a screw-cap bottle designed specifically to meet all state requirements for TCLP extraction and analysis for the metals listed below.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Antimony (Sb) | [7440-36-0] | 0,2 - 200 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 0,5 - 400 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 0,5 - 5000 | mg/L | | |
| | Beryllium (Be) | [7440-41-7] | 0,1 - 60 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 0,5 - 60 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 0,5 - 500 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 0,5 - 1500 | mg/L | | |
| | Mercury (Hg) | [7439-97-6] | 0,05 - 100 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 0,5 - 200 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 0,5 - 100 | mg/L | | |
| | Silver (Ag) | [7440-22-4] | 0,2 - 400 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 0,5 - 300 | mg/L | | |