

**Glass - Setting up sample (SUS) - K-Na, disc 40 mm x 5 mm**

Art. ID XRF-RKU-002  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	SiO <sub>2</sub>		63,65	%		
	Al <sub>2</sub> O <sub>3</sub>		2	%		
	Li <sub>2</sub> O		1	%		
	Na <sub>2</sub> O		12	%		
	K <sub>2</sub> O		9	%		
	MgO	[1309-48-4]	2	%		
	CaO		7	%		
	TiO <sub>2</sub>		0,2	%		
	Cr <sub>2</sub> O <sub>3</sub>		0,2	%		
	MoO <sub>3</sub>		0,5	%		
	Fe <sub>2</sub> O <sub>3</sub>		1	%		
	Co <sub>3</sub> O <sub>4</sub>		0,1	%		
	As <sub>2</sub> O <sub>3</sub>		0,5	%		
	Sb <sub>2</sub> O <sub>3</sub>		0,5	%		
	SO <sub>3</sub>		0,2	%		
	Selenium (Se)	[7782-49-2]	0,05	%		
	MnO		0,1	%		