

**Glass - Setting up sample (SUS) - F-P Glass, disc 40 mm x 5 mm**

Art. ID XRF-RKU-008  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	SiO <sub>2</sub>		40,55	%		
	Al <sub>2</sub> O <sub>3</sub>		6	%		
	Na <sub>2</sub> O		7	%		
	K <sub>2</sub> O		7	%		
	MgO	[1309-48-4]	3	%		
	CaO		5	%		
	BaO		5	%		
	ZrO <sub>2</sub>		0,02	%		
	B <sub>2</sub> O <sub>3</sub>		10	%		
	ZnO		1	%		
	P <sub>2</sub> O <sub>5</sub>		5	%		
	Sb <sub>2</sub> O <sub>3</sub>		0,4	%		
	SO <sub>3</sub>		0,03	%		
	Fluorine (F)	[7782-41-4]	4	%		
	CeO <sub>2</sub>		1	%		
	SrO		5	%		