

Aluminium base, CRM, set of 6 discs each approx 60x30/20mm (DHJ-58A-ZB23013, DHJ-58A-ZB23014, DHJ-58A-ZB23015a, DHJ-58A-ZB23016a, DHJ-58A-ZB23017, DHJ-58A-ZB23018)

Art. ID DHJ-58A-ZB23013-ZB3018
Unit discs (6)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
DHJ-58A-ZB23013: Alumin Beryllium (Be) ium base, CRM, disc 60x 30mm	Beryllium (Be)	[7440-41-7]	0,0061	%		
DHJ-58A-ZB23013: Alumin Calcium (Ca) ium base, CRM, disc 60x 30mm	Calcium (Ca)	[7440-70-2]	0,0027	%		
DHJ-58A-ZB23013: Alumin Cadmium (Cd) ium base, CRM, disc 60x 30mm	Cadmium (Cd)	[7440-43-9]	0,0013	%		
DHJ-58A-ZB23013: Alumin Chromium (Cr) ium base, CRM, disc 60x 30mm	Chromium (Cr)	[7440-47-3]	0,0051	%		
DHJ-58A-ZB23013: Alumin Copper (Cu) ium base, CRM, disc 60x 30mm	Copper (Cu)	[7440-50-8]	0,01	%		
DHJ-58A-ZB23013: Alumin Iron (Fe) ium base, CRM, disc 60x 30mm	Iron (Fe)	[7439-89-6]	0,026	%		
DHJ-58A-ZB23013: Alumin Gallium (Ga) ium base, CRM, disc 60x 30mm	Gallium (Ga)	[7440-55-3]	0,0054	%		
DHJ-58A-ZB23013: Alumin Magnesium (Mg) ium base, CRM, disc 60x 30mm	Magnesium (Mg)	[7439-95-4]	0,196	%		
DHJ-58A-ZB23013: Alumin Manganese (Mn) ium base, CRM, disc 60x 30mm	Manganese (Mn)	[7439-96-5]	0,085	%		
DHJ-58A-ZB23013: Alumin Nickel (Ni) ium base, CRM, disc 60x 30mm	Nickel (Ni)	[7440-02-0]	0,0053	%		

DHJ-58A-ZB23013: Alumin Phosphorus (P)	[7723-14-0]	0,0018	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Lead (Pb)	[7439-92-1]	0,0053	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Silicon (Si)	[7440-21-3]	5,56	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Tin (Sn)	[7440-31-5]	0,0058	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Strontium (Sr)	[7440-24-6]	0,0064	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Titanium (Ti)	[7440-32-6]	0,024	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Vanadium (V)	[7440-62-2]	0,0049	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Zinc (Zn)	[7440-66-6]	0,0064	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23013: Alumin Zirconium (Zr)	[7440-67-7]	0,0047	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23014: Alumin Boron (B)	[7440-42-8]	0,0006	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23014: Alumin Beryllium (Be)	[7440-41-7]	0,01	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23014: Alumin Calcium (Ca)	[7440-70-2]	0,0006	%
ium base, CRM, disc 60x 30mm			
DHJ-58A-ZB23014: Alumin Cadmium (Cd)	[7440-43-9]	0,0061	%
ium base, CRM, disc 60x			

30mm

DHJ-58A-ZB23014: Alumin Chromium (Cr) [7440-47-3] 0,019 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Copper (Cu) [7440-50-8] 0,051 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Iron (Fe) [7439-89-6] 0,214 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Gallium (Ga) [7440-55-3] 0,059 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Magnesium (Mg) [7439-95-4] 0,288 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Manganese (Mn) [7439-96-5] 0,121 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Nickel (Ni) [7440-02-0] 0,045 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Phosphorus (P) [7723-14-0] 0,0062 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Lead (Pb) [7439-92-1] 0,044 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Silicon (Si) [7440-21-3] 7,22 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Tin (Sn) [7440-31-5] 0,084 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Strontium (Sr) [7440-24-6] 0,016 %
ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23014: Alumin Titanium (Ti) [7440-32-6] 0,043 %

ium base, CRM, disc 60x			
30mm			
DHJ-58A-ZB23014: Alumin Vanadium (V)	[7440-62-2]	0,032	%
ium base, CRM, disc 60x			
30mm			
DHJ-58A-ZB23014: Alumin Zinc (Zn)	[7440-66-6]	0,09	%
ium base, CRM, disc 60x			
30mm			
DHJ-58A-ZB23014: Alumin Zirconium (Zr)	[7440-67-7]	0,006	%
ium base, CRM, disc 60x			
30mm			
DHJ-58A-ZB23015a: Alumi Chromium (Cr)	[7440-47-3]	0,048	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Copper (Cu)	[7440-50-8]	0,05	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Iron (Fe)	[7439-89-6]	0,14	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Magnesium (Mg)	[7439-95-4]	0,29	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Manganese (Mn)	[7439-96-5]	0,05	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Nickel (Ni)	[7440-02-0]	0,052	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Lead (Pb)	[7439-92-1]	0,051	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Silicon (Si)	[7440-21-3]	6,85	%
nium base, CRM, disc 50			
x30mm			
DHJ-58A-ZB23015a: Alumi Tin (Sn)	[7440-31-5]	0,049	%
nium base, CRM, disc 50			
x30mm			

DHJ-58A-ZB23015a: Alumi Strontium (Sr) nium base, CRM, disc 50 x30mm	[7440-24-6]	0,019	%
DHJ-58A-ZB23015a: Alumi Titanium (Ti) nium base, CRM, disc 50 x30mm	[7440-32-6]	0,055	%
DHJ-58A-ZB23015a: Alumi Zinc (Zn) nium base, CRM, disc 50 x30mm	[7440-66-6]	0,052	%
DHJ-58A-ZB23016a: Alumi Chromium (Cr) nium base, CRM, disc 50 x30mm	[7440-47-3]	0,022	%
DHJ-58A-ZB23016a: Alumi Copper (Cu) nium base, CRM, disc 50 x30mm	[7440-50-8]	0,052	%
DHJ-58A-ZB23016a: Alumi Iron (Fe) nium base, CRM, disc 50 x30mm	[7439-89-6]	0,117	%
DHJ-58A-ZB23016a: Alumi Magnesium (Mg) nium base, CRM, disc 50 x30mm	[7439-95-4]	0,379	%
DHJ-58A-ZB23016a: Alumi Manganese (Mn) nium base, CRM, disc 50 x30mm	[7439-96-5]	0,053	%
DHJ-58A-ZB23016a: Alumi Nickel (Ni) nium base, CRM, disc 50 x30mm	[7440-02-0]	0,031	%
DHJ-58A-ZB23016a: Alumi Lead (Pb) nium base, CRM, disc 50 x30mm	[7439-92-1]	0,031	%
DHJ-58A-ZB23016a: Alumi Silicon (Si) nium base, CRM, disc 50 x30mm	[7440-21-3]	6,99	%
DHJ-58A-ZB23016a: Alumi Tin (Sn) nium base, CRM, disc 50 x30mm	[7440-31-5]	0,029	%
DHJ-58A-ZB23016a: Alumi Strontium (Sr) nium base, CRM, disc 50	[7440-24-6]	0,016	%

x30mm				
DHJ-58A-ZB23016a: Alumi Titanium (Ti)	[7440-32-6]	0,133		%
nium base, CRM, disc 50				
x30mm				
DHJ-58A-ZB23016a: Alumi Zinc (Zn)	[7440-66-6]	0,024		%
nium base, CRM, disc 50				
x30mm				
DHJ-58A-ZB23017: Alumin Beryllium (Be)	[7440-41-7]	0,12		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Calcium (Ca)	[7440-70-2]	0,017		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Cadmium (Cd)	[7440-43-9]	0,058		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Chromium (Cr)	[7440-47-3]	0,096		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Copper (Cu)	[7440-50-8]	0,191		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Iron (Fe)	[7439-89-6]	0,373		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Gallium (Ga)	[7440-55-3]	0,096		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Magnesium (Mg)	[7439-95-4]	0,497		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Manganese (Mn)	[7439-96-5]	0,207		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Nickel (Ni)	[7440-02-0]	0,067		%
ium base, CRM, disc 60x				
30mm				
DHJ-58A-ZB23017: Alumin Phosphorus (P)	[7723-14-0]	0,0033		%

ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Lead (Pb)	[7439-92-1]	0,101		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Silicon (Si)	[7440-21-3]	6,26		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Tin (Sn)	[7440-31-5]	0,036		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Strontium (Sr)	[7440-24-6]	0,034		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Titanium (Ti)	[7440-32-6]	0,153		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Vanadium (V)	[7440-62-2]	0,175		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Zinc (Zn)	[7440-66-6]	0,148		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23017: Alumin Zirconium (Zr)	[7440-67-7]	0,072		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23018: Alumin Beryllium (Be)	[7440-41-7]	0,21		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23018: Alumin Calcium (Ca)	[7440-70-2]	0,0053		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23018: Alumin Cadmium (Cd)	[7440-43-9]	0,12		%
ium base, CRM, disc 60x 30mm				
DHJ-58A-ZB23018: Alumin Chromium (Cr)	[7440-47-3]	0,108		%
ium base, CRM, disc 60x 30mm				

DHJ-58A-ZB23018: Alumin Copper (Cu) ium base, CRM, disc 60x 30mm	[7440-50-8]	0,275	%
DHJ-58A-ZB23018: Alumin Iron (Fe) ium base, CRM, disc 60x 30mm	[7439-89-6]	0,483	%
DHJ-58A-ZB23018: Alumin Gallium (Ga) ium base, CRM, disc 60x 30mm	[7440-55-3]	0,2	%
DHJ-58A-ZB23018: Alumin Magnesium (Mg) ium base, CRM, disc 60x 30mm	[7439-95-4]	0,578	%
DHJ-58A-ZB23018: Alumin Manganese (Mn) ium base, CRM, disc 60x 30mm	[7439-96-5]	0,307	%
DHJ-58A-ZB23018: Alumin Nickel (Ni) ium base, CRM, disc 60x 30mm	[7440-02-0]	0,091	%
DHJ-58A-ZB23018: Alumin Phosphorus (P) ium base, CRM, disc 60x 30mm	[7723-14-0]	0,0035	%
DHJ-58A-ZB23018: Alumin Lead (Pb) ium base, CRM, disc 60x 30mm	[7439-92-1]	0,138	%
DHJ-58A-ZB23018: Alumin Silicon (Si) ium base, CRM, disc 60x 30mm	[7440-21-3]	7,82	%
DHJ-58A-ZB23018: Alumin Tin (Sn) ium base, CRM, disc 60x 30mm	[7440-31-5]	0,057	%
DHJ-58A-ZB23018: Alumin Strontium (Sr) ium base, CRM, disc 60x 30mm	[7440-24-6]	0,013	%
DHJ-58A-ZB23018: Alumin Titanium (Ti) ium base, CRM, disc 60x 30mm	[7440-32-6]	0,184	%
DHJ-58A-ZB23018: Alumin Vanadium (V) ium base, CRM, disc 60x	[7440-62-2]	0,157	%

30mm

DHJ-58A-ZB23018: Alumin Zinc (Zn)	[7440-66-6]	0,338	%
-----------------------------------	-------------	-------	---

ium base, CRM, disc 60x

30mm

DHJ-58A-ZB23018: Alumin Zirconium (Zr)	[7440-67-7]	0,065	%
--	-------------	-------	---

ium base, CRM, disc 60x

30mm