

Sodium Aminosalicylate dihydrate for equipment qualification

Art. ID	EPY0001816
Unit	EP unit
Deliverydetails	No Dangerous Good /not restricted

Description

Monograph: 20232, 20512, 20532 /// Sodium aminosalicylate dihydrate for equipment qualification CRS batch 1 can only be used to perform Ph. Eur. test 2.5.32. Micro determination of water when proceeding by direct or liquid sample introduction. For micro determination of water with oven, use Amoxicillin trihydrate for performance verification CRS (Y0001521).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Sodium 4-aminosalicylat e dihydrate	Sodium aminosalicylate dihydrate	[6018-19-5]				