

National Bureau of Standards

Certificate of Analyses

Standard Sample 100B Manganese Steel

ANALYST	C	Mn	P		S			Si	Cu	Ni	Cr	V	Mo	N
	Direct combustion	Persulfate-Arsenite	Gravimetric (weighed as $Mg_2P_2O_7$ after removal of arsenic)	Alkali-Molybdate ^a	Gravimetric (direct oxidation and precipitation after reduction of iron)	Combustion iodate titration	Evolution with HCl (1+1) ZnS-Iodine (theoretical sulfur titer) ^b	Perchloric acid dehydration		Weighted as nickel dimethylglyoxime	$FeSO_4$ -KMnO ₄ titration		Colorimetric	Distillation-titration
1.....	0.395	1.89	0.021	0.025	0.030	0.028	0.030	0.213	0.064	0.029	0.063	0.003	0.238	0.004
2.....	.394	1.91	.021	.021	.028	.028		.209	.068	.028	.063	.003	.235	.004
3.....	.407	1.89		.024	.030	.029	.028	^m .21	.063	0.030	.062	0.003	.239	.005
4.....	.401	1.90		.023	.027	.027	.028	^f .207	.063	.029	.064	0.003	.242	^k .005
5.....	.392	1.88	.023	0.023	.029		.027	^m .219	.060	0.030	.061	0.003	.233	
6.....	.392	1.88	.023	0.022	.030	0.030		^m .217	.058	0.029	.058	0.003	.229	
7.....	.397	1.87		^w .024			.026	.196	^x .076	0.034	.068	^k .001	.243	
Averages.....	0.397	1.89	0.022	0.023	0.029	0.028	0.028	0.210	0.064	0.030	0.063	0.003	0.237	0.004
General Average.....	0.397	1.89	0.023		0.028			0.210	0.064	0.030	0.063	0.003	0.237	0.004

^a Precipitated at 40°C, washed with a 1-percent solution of KNO₃, and titrated with alkali standardized by the use of acid potassium phthalate and the ratio 23NaOH:1P.

^b Value obtained by standardizing the titrating solution by means of sodium oxalate through KMnO₄ and Na₂S₂O₄ and use of the ratio 2I:1S.

^c Potentiometric titration.

^d Molybdenum-blue photometric method. See J. Research NBS 26, 405 (1941), RP1386.

^e 1-g sample burned in oxygen at 1,425°C and sulfur dioxide absorbed in starch-iodide solution. Iodine liberated from iodide by titration, during the combustion, with standard KIO₃ solution. Titer based on 93 percent of the theoretical factor.

^f Double dehydration with intervening filtration.

^g Diethyldithiocarbamate photometric method. See J. Research NBS 47, 380 (1951), RP2265.

^h Chromium separated from the bulk of the iron in a 10-g sample by hydrolytic precipitation with NaHCO₃, oxidized with persulfate, and titrated potentiometrically with ferrous ammonium sulfate.

ⁱ Vanadium separated as in (h), oxidized with HNO₃, and titrated potentiometrically with ferrous ammonium sulfate.

^j Sulfuric acid digestion for 4 hr of a 0.5-g sample. See J. Research NBS 43 201 (1949) RP2021.

^k Photometric method.

^l Titrating solution standardized with a standard steel.

^m Sulfuric acid dehydration.

ⁿ Copper-ammonia-complex photometric method.

^o Dimethylglyoxime photometric method.

^p Diphenylcarbazide photometric method.

^q NaHCO₃-FeSO₄-(NH₄)₂S₂O₈-KMnO₄ method.

^r H₂S-CuS-CuO.

^s H₂S-MoS₃-MoO₃.

^t CuS precipitated with Na₂S₂O₄, and copper determined by electrolysis.

^u Dimethylglyoxime precipitate ignited to NiO.

^v Vanadium precipitated with cupferron and determined by the FeSO₄-(NH₄)₂S₂O₈-KMnO₄ method.

^w Molybdenum-blue photometric method.

^x Neocuproine photometric method.

^y Perchloric acid-photometric method.

List of Analysts

1. Ferrous Laboratory, National Bureau of Standards. J. I. Shultz, in charge. Analysis by E. June Maienthal, T. W. Freeman, and E. J. Meros.
2. H. J. Wolthorn, United States Steel Corporation, Fairless Works, Fairless Hills, Pa.
3. J. Henderson, United States Steel Corporation, Edgar Thomson Works, Braddock, Pa.

4. D. J. Hallisey, Jones and Laughlin Steel Corporation, Aliquippa Works, Aliquippa, Pa.
5. Charles C. Kawin Company, Buffalo, N.Y.
6. Charles C. Kawin Company, Chicago, Ill.
7. J. T. Crane, Page Steel and Wire Division, American Chain and Cable Company, Inc., Monessen, Pa.

The steel for the preparation of this standard was furnished by the United States Steel Corporation.
WASHINGTON, D.C., August 18, 1959.

A. V. ASTIN, *Director*