



Certificate of Analysis

Standard Reference Material 152A

Basic Open-Hearth Steel 0.5% Carbon

(Tin-Bearing)

ANALYST	C	Mn	P	S	Si	Cu	Ni	Cr	V	Mo	Sn	
	Direct combustion	Persulfate-Arsenite	Photometric	Combustion Iodate titration	Perchloric acid dehydration	Photometric	Weighed as nickel dimethyl- glyoxime	FeSO ₄ -KMnO ₄ titration		Photometric		
1.....	0.484	0.714	^a 0.012	^b 0.030	^c 0.202	^d 0.025	^e 0.057	^f 0.047	^g 0.001	0.035	^h 0.035
2.....	.486	.718	.012	{ ⁱ 0.030 ^k 0.030}	^c .205	^l .022	.055	^m .046	ⁿ <.001	.039	^h .031
3.....	.488	.716	.012	.031	^c .204	^o .021	^e .058	^p .049	^q <.001	.037	.031
4.....	^r .484	.720	.012	.030	.203	^o .024	.052	.041	ⁿ .001	.036	^h .033
5.....	.491	^s .72	.012	.028	.198	^t .023	^e .054037
6.....	.485	.714	.013	.031	.202	^t .021	^e .057	^p .048	ⁿ .003	.033	.031
Average.....	0.486	0.717	0.012	0.030	0.202	0.023	0.056	0.046	0.001	0.036	0.032

^a Molybdenum-blue photometric method. See J. Res. NBS 26, 405 (1941) RP 1386.

^b 1-g sample burned in oxygen at 1425 °C. and sulfur dioxide absorbed in starch-iodide solution. Iodine is liberated from iodide by titration, during the combustion, with standard KIO₃ solution. Titer is based on 93 percent of the theoretical factor.

^c Double dehydration with intervening filtration.

^d Diethyldithiocarbamate photometric method. See J. Res. NBS 47, 380 (1951) RP2265.

^e Dimethylglyoxime photometric method.

^f Chromium separated from the bulk of the iron in a 10-g sample by hydrolytic precipitation with NaHCO₃, oxidized with persulfate, and titrated potentiometrically with ferrous ammonium sulfate.

^g Vanadium separated as in (f), oxidized with HNO₃, and titrated potentiometrically with ferrous ammonium sulfate.

^h Sulfide-iodine method. See BS. J. Res. 8, 309 (1932) RP415.

ⁱ Titrating solution standardized by use of a standard steel.

^j Alkali-molybdate method.

^k Gravimetric method. Sulfate precipitated with BaCl₂, ignited to BaSO₄ and weighed.

^l H₂S-electrolytic method.

^m Sodium bicarbonate hydrolysis-persulfate oxidation.

ⁿ Sodium bicarbonate hydrolysis-FeSO₄-(NH₄)₂S₂O₈-KMnO₄.

^o Copper-ammonia complex photometric method.

^p Diphenylcarbazide photometric method.

^q Spectrographic method.

^r Gasometric method.

^s KIO₃ photometric method.

^t Diethyldithiocarbamate photometric method.

List of Analysts

1. B. B. Bendigo, E. R. Deardorff and J. I. Shultz, Division of Analytical Chemistry, National Bureau of Standards.
2. A. C. Nickel, United States Steel Corp., Gary Works, Gary, Ind.
3. G. K. Stewart, United States Steel Corp., Geneva Works, Provo, Utah.

4. C. A. Spellman, Joslyn Stainless Steels, Fort Wayne, Ind.
5. S. M. Dibble, Crucible Steel Company of America, Syracuse, New York.
6. H. W. Huston, A. M. Byers Co., Ambridge, Pa.

The steel for the preparation of this standard was furnished by the United States Steel Corporation.

October 11, 1965

W. Wayne Meinke, Chief
Office of Standard Reference Materials