



# CERTIFIED REFERENCE MATERIAL BCR-326

## CERTIFICATE OF ANALYSIS

UNALLOYED ZINC			
	Mass Fraction		Number of accepted sets of data p
	Certified value <sup>1)</sup> [µg/g]	Uncertainty [µg/g]	
Cd	203.0	2.0 <sup>2)</sup>	13
Cu	104.8	2.7 <sup>3)</sup>	13
Fe	264.8	2.1 <sup>2)</sup>	11
Pb	307.0	1.6 <sup>2)</sup>	12
<p><sup>1)</sup> Unweighted mean value of p sets of data, each set being obtained in a different laboratory and/or with a different method of determination. The certified values are traceable to the International System of Units (SI).</p> <p><sup>2)</sup> Half-width of the 95 % confidence interval of the certified value as defined in 1).</p> <p><sup>3)</sup> Half-width of the 95 % tolerance interval.</p>			

This certificate is valid for five years after purchase.

Sales date:

The minimum amount of sample to be used is 1 g. It is recommended not to use the center of the sample (10 mm diameter).

### NOTE

This material has been certified by BCR (Community Bureau of Reference, the former reference materials programme of the European Commission). The certificate has been revised under the responsibility of IRMM.

Brussels, November 1989

Latest revision: August 2015

Signed: \_\_\_\_\_

Prof. Dr. Hendrik Emons  
European Commission  
Joint Research Centre  
Institute for Reference Materials and Measurements  
Retieseweg 111  
B-2440 Geel, Belgium

Additional Material Information	
	Mass Fraction
	Value <sup>1</sup> [µg/g]
Al	< 0.7
1) Upper limit of the mean supported by 4 sets of data, each set being obtained in a different laboratory and/or with a different method of determination.	

## DESCRIPTION OF THE SAMPLE

The samples are discs with 80 mm diameter and 20 mm thickness.

Each sample has two different numbers:

- a five digit code, marked on the metal by the producer of the materials, defining the exact position of the sample in the original batch of rods;
- the "normal" individual identification (consecutive numbering of samples within a CRM), marked on the label of the samples.

	Individual identification (on label)	Five digit code (on metal sample) *
BCR-326	001 to 060	16101 to 16160
	061 to 120	16201 to 16260
	121 to 180	16301 to 16360
	181 to 240	16401 to 16460

- \* First digit : 1 for unalloyed zinc  
 Second digit : 6 for BCR-326  
 Third digit : rod number within one CRM  
 Last two digits : disk number within one rod (01 = bottom, 60 = top)

## ANALYTICAL METHOD USED FOR CERTIFICATION

- Differential pulse voltammetry
- Electrothermal atomic absorption spectrometry
- Flame atomic absorption spectrometry
- Inductively coupled plasma emission spectrometry
- Polarography
- Spectrophotometry

## PARTICIPANTS

- Budelco B.V., Budel-Dorplein (NL)
- Bundesanstalt für Materialforschung und –prüfung, Berlin (DE)
- Energieonderzoek Centrum Nederland, Petten (NL)
- Ever Ready Ltd, Durham (GB)
- Instituut voor Nucleaire Wetenschappen, Gent (BE)
- Laboratoire National d'Essais, Paris (FR)
- Metaleurop Harz-Metall GmbH, Goslar (DE)
- Metaleurop, Noyelles-Godault (FR)
- Metaleurop Weiser Blei GmbH, Nordenham (DE)
- Minemet Recherche, Trappes (FR)
- Pertusola Sud, Crotone (IT)
- Prof. G. Kraft, Kronberg (DE)
- Ruhr-Zink GmbH, Datteln (DE)
- Vieille-Montagne, Balen (BE)
- Vieille-Montagne, Overpelt (BE)
- Vieille-Montagne France, Aubry-les-Douai (FR)

## **SAFETY INFORMATION**

Not applicable.

## **INSTRUCTIONS FOR USE**

The material is mainly intended for calibration in emission spectrometry with solid samples; the usual mechanical cleaning should be applied prior to the measurement (the CRM and the user's samples should be treated in the same way). It is recommended not to use the centre of the samples (approx. 10 mm diameter).

## **STORAGE**

The CRMs should be stored in clean and dry conditions at ambient temperature.

The European Commission cannot be held responsible for changes that happen during storage of the material at the customer's premises, especially of opened samples.

## **LEGAL NOTICE**

Neither IRMM, its subsidiaries, its contractors nor any person acting on their behalf,

(a) make any warranty or representation, express or implied that the use of any information, material, apparatus, method or process disclosed in this document does not infringe any privately owned intellectual property rights; or

(b) assume any liability with respect to, or for damages resulting from, the use of any information, material, apparatus, method or process disclosed in this document save for loss or damage arising solely and directly from the negligence of IRMM or any of its subsidiaries.

## **NOTE**

A technical report on the production of BCR-326 is available on the internet (<http://www.irmm.jrc.be>). A paper copy can be obtained from IRMM on request.

---

European Commission – Joint Research Centre  
Institute for Reference Materials and Measurements (IRMM)  
Retieseweg 111, B - 2440 Geel (Belgium)  
Telephone: +32-(0)14-571.722 - Telefax: +32-(0)14-590.406