



CERTIFIED REFERENCE MATERIAL BCR[®] – 420

CERTIFICATE OF ANALYSIS

WASTE MINERAL OIL				
Congener number	IUPAC name	Mass Fraction based on dry mass		Number of accepted sets of results p
		Certified value [mg/kg] ¹⁾	Uncertainty [mg/kg] ²⁾	
28	2,4,4'-Trichlorobiphenyl	0.61	0.06	10
101	2,2',4,5,5'-Pentachlorobiphenyl	1.45	0.18	5
118	2,3',4,4',5-Pentachlorobiphenyl	1.69	0.14	6
153	2,2',4,4',5,5'-Hexachlorobiphenyl	0.92	0.06	10
180	2,2',3,4,4',5,5'-Heptachlorobiphenyl	0.195	0.017	10

¹⁾ This value is the unweighted mean of the means of p accepted sets of results. The certified values are traceable to determinations by GC-based methods.
²⁾ The uncertainty is taken as the half-width of the 95 % confidence interval of the mean given in ¹⁾.

This certificate is valid for one year after purchase.

Sales date:

The minimum amount of sample to be used is 200 mg.

NOTE

This material has been certified by BCR (Community Bureau of Reference, the former reference materials programme of the European Commission). The certificate has been revised under the responsibility of IRMM.

Brussels, January 1991

Revised: April 2007

Signed: _____

Prof. Dr. Hendrik Emons
Unit for Reference Materials
EC-JRC-IRMM
Retieseweg 111
2440 Geel, Belgium

DESCRIPTION OF THE SAMPLE

The sample consists of approximately 7.5 g of waste mineral oil in brown glass ampoules sealed under argon. Additional information on the presence of chlorobiphenyls other than those tabulated is given in the report.

ANALYTICAL METHOD USED FOR CERTIFICATION

Calibration was done with solutions of BCR-CRMs (291, 293-298). The samples were dissolved in hexane and clean-up was carried out by sulfuric acid and/or column chromatography (on alumina, silica gel or florisil). Capillary gas chromatography with electron capture detection was performed using different injection systems, different columns and different temperature programmes.

PARTICIPANTS

- Department of Agriculture and Fisheries for Scotland (DAFS), Aberdeen (GB)
- Direcção General da Qualidade do Ambiente, Lisbon (PT)
- European Commission, Joint Research Centre, Environment Institute, Ispra (IT)
- Institute for Environmental Studies, Free University of Amsterdam, Amsterdam (NL)
- Instituto Químico de Sarria, Barcelona (ES)
- Institute of Marine Research, Bergen (NO)
- Instituto Hidrográfico, Lisbon (PT)
- Laboratoire Municipal et Régional, Rouen (FR)
- National Swedish Environment Protection Board, Solna (SE)
- Netherlands Institute for Fishery Invest. (RIVO), IJmuiden (NL)
- Studiecentrum voor Kernenergie SCK/CEN, Mol (BE)

SAFETY INFORMATION

Cancer suspect agent. The material must be handled with great care, especially avoiding skin contamination, ingestion or inhalation.

INSTRUCTIONS FOR USE

For analysis the sample should be taken as it is. The water content is typically below 1 % by mass. In case the user would need to correct for the water content a titration is recommended (e.g. Karl Fischer).

STORAGE

After opening, the content should be transferred to a clean glass container which has a ground glass stopper, then stored at 5 °C in the dark.

However, the European Commission cannot be held responsible for changes that happen during storage of the material at the customer's premises, especially of opened samples.

LEGAL NOTICE

Neither IRMM, its subsidiaries, its contractors nor any person acting on their behalf.

(a) make any warranty or representation, express or implied that the use of any information, material, apparatus, method or process disclosed in this document does not infringe any privately owned intellectual property rights; or

(b) assume any liability with respect to, or for damages resulting from, the use of any information, material, apparatus, method or process disclosed in this document save for loss or damage arising solely and directly from the negligence of IRMM or any of its subsidiaries.

NOTE

A technical report on the production of BCR-420 is available on the internet (<http://www.irmm.jrc.be>). A paper copy can be obtained from IRMM on request.