

**Out of Stock - Item is not available at this time - Gasoline (Reformulated) - Sulfur**

Art. ID                      NIST-2299  
Unit                         5 x 20 mL  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1203 II E2

Description

This Standard Reference Material (SRM®) is a commercial reformulated unleaded motor gasoline intended for use in the evaluation of methods and the calibration of instruments used in the determination of total sulfur in gasoline or materials of similar matrix. A unit of NIST-2299 consists of 5 amber ampoules, each containing approximately 20 mL of gasoline. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	13,6 ± 1,5	µg/g		