

## **Titanium Alloy**

Art. ID                      NIST-2433  
Unit                         50 g  
Deliverydetails            No Dangerous Good /not restricted

### **Description**

This Standard Reference Material (SRM®) is in the form of chips sized between 0.50 and 1.18 mm sieve openings (35 and 16 mesh). It is intended for use primarily in chemical methods analysis. Certified values /// Sample value(s)  
- please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	7,63	%		
	Iron (Fe)	[7439-89-6]	0,063	%		
	Molybdenum (Mo)	[7439-98-7]	0,99	%		
	Vanadium (V)	[7440-62-2]	0,98	%		