

Lead paint film (white/blank) - Lead (Nominal <0.001 mg/cm²)

Art. ID	NIST-2570
Unit	1 film
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is intended as a blank for NIST-2571 to NIST-2576, which are used for checking the calibration of portable, hand-held, x-ray fluorescence analyzers when testing for lead in paint coatings on interior and exterior building surfaces. A unit of NIST-2570 consists of a white polyester sheet, approximately 7.6 cm wide, 10.2 cm long, and 0.2 mm thick, coated with a single, lacquer layer, approximately 0.04 mm thick. Each unit is over-coated with a clear, thin, plastic laminate to protect the surface from abrasion. NIST-2570 is one of a set of six paint films (NIST-2570 to NIST-2575) available as NIST-2579a. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Color code: White (Blank)	Lead (Pb)	[7439-92-1]	<0,001	mg/cm ²		