

**Out of Stock - Product is not available at this time - Water saturated octanol**

Art. ID	NIST-2890
Unit	5 x 2 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Standard Reference Material (SRM®) 2890 is intended primarily for use in validating methods and calibrating instruments for the measurement of trace concentration levels of water. NIST-2890 is a solution of water saturated 1-octanol which is certified for its water content. A unit of NIST-2890 consists of five 5 mL ampoules, each containing approximately 2 mL of water saturated 1-octanol. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Water content		47,3 ± 1,0	mg/g		
	Water content		39,24 ± 0,85	mg/mL		