

**Out of Stock - Item is not available at this time - Human cardiac Troponin complex**

Art. ID                      NIST-2921  
Unit                         5 x 115 uL  
Deliverydetails            Dry ice shipment /not restricted

Description

This Standard Reference Material (SRM®) is primarily intended for use in calibrating clinical procedures and devices for the determination of cardiac troponin I (cTnI) in human serum. It can also be used for value-assignment of calibrators and control materials. A unit of NIST-2921 consists of five vials, each containing approximately 115 µL of a dilute solution of human cardiac troponin complex. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Certified cTnI Concentration		31,2 ± 1,4	mg/L		