

Arsenic (As) Standard Solution (nominal 10 mg/g)

Art. ID	NIST-3103a
Unit	5 x 10 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 II E2

Description

Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
As in 15% HNO ₃ , 10 mg/g	Arsenic (As)	[7440-38-2]	10	mg/g		