

Beryllium (Be) Standard Solution (nominal 10 mg/g)

Art. ID NIST-3105a
Unit 5 x 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 II E2

Description

Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Be in 10% HNO ₃ , Lot No. 090514	Beryllium (Be)	[7440-41-7]	9,960 ± 0,034	mg/g		