

Praseodymium (Pr) Standard Solution

Art. ID	NIST-3142a
Unit	5 x 10 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 II E2

Description

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of praseodymium. One unit of the SRM consists of five 10 mL sealed borosilicate glass ampoules of an acidified aqueous solution prepared gravimetrically to contain a known mass fraction of praseodymium. The solution contains nitric acid at a mass fraction of approximately 17 %. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Praseodymium (Pr)	[7440-10-0]	10,012 ± 0,023	mg/g		