

## **Sodium (Na) Standard Solution**

Art. ID                      NIST-3152a  
Unit                        50 mL  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 II E2

### **Description**

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of sodium. One unit of NIST-3152a consists of 50 mL of an aqueous solution in a high-density polyethylene bottle sealed in an aluminized bag. The solution was prepared gravimetrically to contain a known mass fraction of sodium in nitric acid at a mass fraction of 5 %. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Lot No. 200413	Sodium (Na)	[7440-23-5]	10,005 ± 0,019	mg/g		