

## **Human Liver (Radioactivity)**

Art. ID	NIST-4352
Unit	45 g
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) consists of approximately 45 grams of dried, cryogenically ground human liver (2), vacuum packed in a 125 ml glass bottle. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	238Pu		5,5 x 10 <sup>-5</sup>	Bq/g		
	239Pu + 240Pu		2,06 x 10 <sup>-3</sup>	Bq/g		
	241Am		1,5 x 10 <sup>-4</sup>	Bq/g		