

## Lake sediment

Art. ID	IAEA-SL-2
Unit	250 g
Deliverydetails	No Dangerous Good /not restricted

### Description

This sample is intended to be used as a reference material for the measurement of 40K and 137Cs in lake sediment samples. It can also be used as a quality control material for the assessment of a laboratory's analytical work, for the validation of analytical methods and for quality assurance within a laboratory.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Recommended Value (Bas: 137Cs d on dry weight)			2,4	Bq/kg		
Recommended Value (Bas: 40K d on dry weight)			240	Bq/kg		