

Ginsenoside - Mix 7 je 500 µg/ml in Ethanol

Art. ID LM24-CUS-10081
Unit 1 mL
Deliverydetails De Minimis UN1170 III E1

Description

Haltbarkeit ungeöffnet 12 Monate

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|-------|--------|--------|
| | Ginsenoside Rg1 | [22427-39-0] | 500 | µg/mL | | |
| | Ginsenoside Re | [52286-59-6] | 500 | µg/mL | | |
| | Ginsenoside Rf | [52286-58-5] | 500 | µg/mL | | |
| | Ginsenoside Rb1 | [41753-43-9] | 500 | µg/mL | | |
| | Ginsenoside Rd | [52705-93-8] | 500 | µg/mL | | |
| | Ginsenoside Rc | [11021-14-0] | 500 | µg/mL | | |
| | Ginsenoside Rb2 | [11021-13-9] | 500 | µg/mL | | |