

## **Commercial Bronze C (block)**

Art. ID                      NIST-C1117  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is designed for use in optical emission and x-ray fluorescence and spectrochemical methods of analysis. Size: 1 1/4 inch in diameter, 3/4 inch thick. Certified values /// Sample value(s) - please ask for current certificate.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Copper (Cu)       | [7440-50-8] | 93,01               | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 6,87                | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,014               | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,02                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,02                | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,069               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,002               | %    |        |        |