

## **Bullet Lead**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | NIST-C2416                        |
| Unit            | disc                              |
| Deliverydetails | No Dangerous Good /not restricted |

## **Description**

This Standard Reference Material (SRM®) is in the form of a disk 50 mm (2 in) in diameter and 16 mm (5/8 in) thick. It is primarily intended for use in X-ray and optical emission spectrometric methods of analysis. /// Sample value(s) - please ask for current certificate.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Antimony (Sb)     | [7440-36-0] | 0,79                | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,056               | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,1                 | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,065               | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,0015              | %    |        |        |
|                  | Silver (Ag)       | [7440-22-4] | 0,0044              | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,09                | %    |        |        |