

## Lead-Base Alloy

Art. ID	NIST-C2417
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is in the form of a disk 50 mm (2 in) in diameter and 16 mm (5/8 in) thick. It is primarily intended for use in X-ray and optical emission spectrometric methods of analysis. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	0,01	%		
	Arsenic (As)	[7440-38-2]	0,011	%		
	Bismuth (Bi)	[7440-69-9]	0,01	%		
	Copper (Cu)	[7440-50-8]	0,01	%		
	Silver (Ag)	[7440-22-4]	0,01	%		