

1-Acetylpiperazine

Art. ID TRC-A187245-100G
Unit 100 g
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) /not restricted

Description

Application Note: 1-Acetylpiperazine was used to study the physical stabilization or chemical degradation of concentrated solutions of polyaniline emeraldine base containing secondary amine additives. It was also used to study the structure-property relationships of surfaces that resist protein adsorption

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|--------------|---------------------|------|--------|--------|
| | 1-Acetylpiperazine | [13889-98-0] | | | | |