

4-Acetylphenylboronic Acid

Art. ID TRC-A188840-5G

Unit 5 g

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: 4-Acetylphenylboronic acid4-Acetylphenylboronic Acid is used in several metal catalyzed cross-coupling reaction studies.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Acetylphenylboronic a cid	[149104-90-5]				