

Acrolein Diethyl Acetal

Art. ID TRC-A191205-2.5G
Unit 2,5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2374 II

Description

Application Note: Acrolein Diethyl Acetal is used as a reagent in quinoline synthesis by improved Skraup-Doebner-Von Miller reactions.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acrolein diethyl acetal	[3054-95-3]				