

AD-mix-beta (Technical Grade)

Art. ID TRC-A200005-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Application Note: AD-mix-beta is a catalyst employed in the Sharpless Asymmetric Dihydroxylation of (E,E)- or (E,Z)-1,3-dienoates. It is also used in the preparation of β -hydroxy- γ -lactones from β , γ -unsaturated esters.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AD-mix-beta	[148618-32-0]				