

## **Adipic Acid-d4**

Art. ID TRC-A291594-50MG  
Unit 50 mg  
Deliverydetails De Minimis

### Description

Application Note: Hexanedioic-3,3,4,4-d4 Acid is the labeled version of A291590. Adipic Acid is primarily used in the synthesis of nylon. It has been used as a reagent in the solid-state polymerization of nylon analogs.

| Text/Information | Analyte/Parameter | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|------|--------|--------|
|                  | Adipic Acid-d4    | [121311-78-2] |                     |      |        |        |