

L-Alanyl-L-glutamine

Art. ID TRC-A511370-5G

Unit 5 g

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: L-Alanyl-L-glutamine is a stable synthetic glutamine dipeptide. This dipeptide is used as an alternative source of amino acids in cell culture medium.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Alanyl-L-Glutamine	[39537-23-0]				