

4-((4-(Allyloxy)phenyl)sulfonyl)phenol

Art. ID TRC-A569030-250MG
Unit 250 mg
Deliverydetails De Minimis /not restricted

Description

Application Note: 4-((4-(Allyloxy)phenyl)sulfonyl)phenol is an additive to enhance the flame-retardation of unsaturated polyester resin (FR-UPR).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-((4-(Allyloxy)phenyl)sulfonyl)phenol	[97042-18-7]				