

[2-(2-Aminoethoxy)ethoxy]acetic Acid

Art. ID TRC-A609435-500MG

Unit 500 mg

Deliverydetails De Minimis

Description

Application Note: [2-(2-Aminoethoxy)ethoxy]acetic Acid is used in the development of chemically modified peptide nucleic acid (PNA) as a probe for qualitative and quantitative detection of DNA. It is also used in the preparation and studies of the biological activities of kappa agonist CovX-bodies.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	[2-(2-Aminoethoxy)ethoxy]acetic Acid	[134978-97-5]				