

2-Amino-4-methylimidazole

Art. ID TRC-A616270-2.5G
Unit 2,5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Application Note: 2-Amino-4-methylimidazole is used as a reagent in the synthesis of 1H-pyrrolo[2,3-c]pyridine-7-carboxamides as novel allosteric mGluR5 antagonists. Also used as a reagent in the synthesis of imidazo[1,2-a]pyrimidin-5(1H)-ones as beta isoform selective phosphatidylinositol 3-kinase (PI3K) inhibitors.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-4-methylimidazo le	[6653-42-5]				