

5-Amino-6-nitroquinoline

Art. ID TRC-A618555-2.5G
Unit 2,5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: 5-Amino-6-nitroquinoline is a reagent used in the synthesis of benzodiazepine sulfonamide-based bombesin receptor 3 agonists. It is also used in the synthesis of quinoxalinediones.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	5-Amino-6-nitroquinolin e	[35975-00-9]				