

Ammonium Sulfate (2:1)

Art. ID TRC-A633835-25G
Unit 25 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: Ammonium Sulfate (2:1) is used in a salting-out method for the purification of tag-free peptides containing multiple disulfide bond in *Pichia pastoris*.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium sulfate	[7783-20-2]				