

AMP PNP

Art. ID TRC-A634303-2.5MG
Unit 2,5 mg
Deliverydetails De Minimis /not restricted

Description

Application Note: AMP PNP is an analogue of AMP (A281780), which is a naturally occurring nucleotide and a useful ligand determinant that facilitates the binding of APS reductase inhibitors and activates adenosine receptor agonists.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ADENYLYL-IMIDODIPHOS [72957-42-7]					
	TE (AMP-PNP),LI4					