

**pHBuffer Solution pH7.000 (Yellow) ±0.010 @25°C (High Resolution)**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | REA-H70525C                       |
| Unit            | 500 mL                            |
| Deliverydetails | No Dangerous Good /not restricted |

**Description**

Coloured solutions. Tested at 20°C or 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory. Traceable to a NIST Standard. Tolerance ± 0.010 pH