

**pHBuffer Capsules pH10.00 (Blue)  $\pm 0.02$  @25°C**

Art. ID	REA-CP1100
Unit	50pk
Deliverydetails	No Dangerous Good /not restricted

**Description**

PH buffers presented in innovative capsule format developed by Reagecon. Colour coded. Easy to use.  
Preservative free. Extended shelf life. Tested at 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory.  
Traceable to a NIST Standard. Accuracy  $\pm 0.05$