

**pHBuffer Capsules Mixed Pack 10xpH4.01, 20xpH7.00, 10xpH9.00, 10xpH10.00 ±0.02
@25°C**

Art. ID	REA-CPMX47910
Unit	50pk
Deliverydetails	No Dangerous Good /not restricted

Description

PH buffers presented in innovative capsule format developed by Reagecon. Colour coded. Easy to use.
Preservative free. Extended shelf life. Tested at 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory.
Traceable to a NIST Standard. Accuracy ±0.06