

**pHBuffer Solution pH7.00 (Yellow)  $\pm 0.01$  @20°C (RECAL)**

Art. ID	REA-07C60
Unit	6 x 90 mL
Deliverydetails	No Dangerous Good /not restricted

**Description**

Supplied in convenient and economical wide mouth, single-use, 6x90ml container designed for Field and Lab direct calibration. Tested at 20 or 25°C and Certified in Reagecon ISO 17025 Accredited Laboratory. Traceable to a NIST Standard. Accuracy  $\pm 0.02$